



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

Product Name:**Malignant melanoma of skin tissue array**

Catalog Number:

TMA- 803

Total Cases//Cores:

39/81

Storage Conditions:

Cool Dry

Array Panel:


Duplicated cores per case of 36 cases of malignant melanoma and 3 cases of adjacent tissue , single care per case

Array Core Diameter:

1.0 mm

Section Thickness:

5 micrometers**TISSUE CORE PANEL DISPLAY ON SLIDE**

		1	2	3	4	5	6	7	8	9
	A									
	B									
	C									
	D									
	E									
	F									
	G									
	H									
	I									

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	88	F	Skin	Malignant melanoma of sole of foot	-	
A2	88	F	Skin	Epidermis	-	
A3	57	M	Skin	Malignant melanoma of scapula	-	
A4	57	M	Skin	Malignant melanoma of scapula	-	
A5	54	F	Skin	Malignant melanoma of sole of foot	-	
A6	54	F	Skin	Malignant melanoma of sole of foot	-	
A7	63	M	Skin	Malignant melanoma of chest wall	-	
A8	63	M	Skin	Malignant melanoma of chest wall	-	
A9	42	F	Skin	Malignant melanoma of thumb	-	
B1	42	F	Skin	Malignant melanoma of thumb	-	
B2	58	F	Skin	Malignant melanoma of buttock	-	
B3	58	F	Skin	Malignant melanoma of buttock	-	
B4	38	M	Skin	Malignant melanoma of abdominal wall	-	
B5	38	M	Skin	Malignant melanoma of abdominal wall	-	
B6	47	F	Skin	Malignant melanoma of anus	-	
B7	47	F	Skin	Malignant melanoma of anus	-	
B8	62	M	Skin	Malignant melanoma of heel	-	
B9	62	M	Skin	Malignant melanoma of heel	-	
C1	82	M	Skin	Malignant melanoma of face	-	
C2	82	M	Skin	Malignant melanoma of face	-	
C3	52	F	Skin	Malignant melanoma of vulvae	-	
C4	52	F	Skin	Malignant melanoma of vulvae	-	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C5	62	F	Skin	Malignant melanoma of thumb	-	
C6	62	F	Skin	Malignant melanoma of thumb	-	
C7	66	F	Skin	Malignant melanoma of chest wall	-	
C8	66	F	Skin	Malignant melanoma of chest wall	-	
C9	25	M	Skin	Malignant melanoma of heel	-	
D1	25	M	Skin	Malignant melanoma of heel	-	
D2	62	M	Skin	Malignant melanoma of thumb	-	
D3	62	M	Skin	Malignant melanoma of thumb	-	
D4	38	F	Skin	Malignant melanoma of toes	-	
D5	38	F	Skin	Malignant melanoma of toes	-	
D6	26	F	Vulva	Malignant melanoma of vulvae	-	
D7	26	F	Vulva	Malignant melanoma of vulvae	-	
D8	52	M	Skin	Malignant melanoma of heel	-	
D9	52	M	Skin	Malignant melanoma of heel	-	
E1	49	M	Skin	Malignant melanoma of shoulder	-	
E2	49	M	Skin	Malignant melanoma of shoulder	-	
E3	36	M	Skin	Malignant melanoma of shoulder	-	
E4	36	M	Skin	Malignant melanoma of shoulder	-	
E5	71	M	Skin	Malignant melanoma of thigh	-	
E6	71	M	Skin	Malignant melanoma of thigh	-	
E7	46	F	Skin	Malignant melanoma of thigh	-	
E8	46	F	Skin	Malignant melanoma of thigh	-	
E9	35	M	Skin	Malignant melanoma of dorsum of foot	-	
F1	35	M	Skin	Malignant melanoma of dorsum of foot	-	
F2	66	M	Skin	Malignant melanoma of thigh	-	
F3	66	M	Skin	Malignant melanoma of thigh	-	
F4	75	F	Skin	Malignant melanoma of hand	-	
F5	75	F	Skin	Malignant melanoma of hand	-	
F6	81	F	Skin	Malignant melanoma of arm	-	
F7	81	F	Skin	Malignant melanoma of arm	-	
F8	65	M	Skin	Malignant melanoma of epicranium	-	
F9	65	M	Skin	Malignant melanoma of epicranium	-	
G1	45	M	Skin	Malignant melanoma of anus	-	
G2	45	M	Skin	Malignant melanoma of anus	-	
G3	48	M	Skin	Malignant melanoma of thumb	-	
G4	48	M	Skin	Malignant melanoma of thumb	-	
G5	45	F	Vulva	Malignant melanoma of vulvae	-	
G6	45	F	Vulva	Malignant melanoma of vulvae	-	
G7	72	F	Vulva	Malignant melanoma of vulvae	-	
G8	72	F	Vulva	Malignant melanoma of vulvae	-	
G9	57	M	Skin	Malignant melanoma of cheeks	-	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

H1	57	M	Skin	Malignant melanoma of cheeks	-	
H2	53	F	Skin	Epidermis	-	
H3	53	F	Skin	Malignant melanoma of thumb	-	
H4	47	F	Skin	Malignant melanoma of upper arm	-	
H5	47	F	Skin	Malignant melanoma of upper arm	-	
H6	53	M	Skin	Malignant melanoma of sole of foot	-	
H7	53	M	Skin	Malignant melanoma of sole of foot	-	
H8	41	F	Skin	Malignant melanoma of groin	-	
H9	41	F	Skin	Malignant melanoma of groin	-	
I1	61	M	Skin	Malignant melanoma of sole of foot	-	
I2	61	M	Skin	Malignant melanoma of sole of foot	-	
I3	59	F	Skin	Malignant melanoma of upper arm	-	
I4	59	F	Skin	Malignant melanoma of upper arm	-	
I5	51	F	Skin	Fibrous tissue and fatty tissue	-	
I6	51	F	Skin	Fibrous tissue and fatty tissue	-	
I7	75	F	Skin	Skin tissue	-	
I8	48	F	Skin	Skin tissue	-	
I9	16	M	Skin	Intradermal nevus	-	

For Research Use Only; Not For Diagnostic or Therapeutic Use.

Conditions: Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA