


Product Name: Breast carcinoma tissue array
Catalog Number: TMA-552
Total Cases//Cores: 80/80
Storage Conditions: Cool Dry
Array Panel: Single core per case, 72 cases of carcinoma and 8 cases of normal tissue
Array Core Diameter: 1.5 mm
Section Thickness: 5 micrometers

TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8
	A								
	B								
	C								
	D								
	E								
	F								
	G								
	H								
	I								
	J								

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	28	F	Breast	Breast tissue	-	
A2	51	F	Breast	Breast tissue	-	
A3	47	F	Breast	Breast tissue	-	
A4	37	F	Breast	Fibrous tissue	-	
A5	44	F	Breast	Breast tissue	-	
A6	45	F	Breast	Breast tissue	-	
A7	38	F	Breast	Fibrofatty tissue	-	
A8	49	F	Breast	Breast tissue	-	
B1	28	F	Breast	Medullary carcinoma	-	
B2	35	F	Breast	Medullary carcinoma	-	
B3	48	F	Breast	Medullary carcinoma	-	
B4	51	F	Breast	Atypical medullary carcinoma with necrosis	-	
B5	47	F	Breast	Mucinous carcinoma	-	
B6	38	F	Breast	Infiltrating lobular carcinoma	-	
B7	38	F	Breast	Infiltrating lobular carcinoma	-	
B8	59	F	Breast	Infiltrating lobular carcinoma	-	
C1	50	F	Breast	Tubular carcinoma	-	
C2	50	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
C3	46	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
C4	36	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
C5	35	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
C6	42	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
C7	40	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C8	49	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
D1	69	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D2	40	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D3	55	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
D4	38	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
D5	52	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
D6	40	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D7	46	F	Breast	A little nonspecific infiltrating ductal carcinoma	-	
D8	53	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E1	46	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E2	53	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
E3	56	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
E4	64	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
E5	52	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E6	56	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E7	35	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E8	69	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F1	70	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F2	52	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F3	49	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F4	40	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F5	52	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F6	57	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F7	62	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F8	38	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G1	50	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G2	38	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G3	46	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
G4	54	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G5	32	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G6	41	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G7	65	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G8	63	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
H1	45	F	Breast	A little nonspecific infiltrating ductal carcinoma	Moderately	
H2	39	F	Breast	A little nonspecific infiltrating ductal carcinoma	Moderately	
H3	41	F	Breast	A little nonspecific infiltrating ductal carcinoma	Moderately	
H4	51	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
H5	42	F	Breast	A little nonspecific infiltrating ductal carcinoma	Moderately	
H6	30	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
H7	53	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
H8	47	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

I1	51	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
I2	46	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
I3	50	F	Breast	Breast tissue	-	
I4	38	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
I5	56	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
I6	45	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
I7	48	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
I8	95	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
J1	59	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
J2	35	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
J3	50	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
J4	50	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
J5	45	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
J6	61	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
J7	74	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
J8	38	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	

For Research Use Only; Not For Diagnostic or Therapeutic Use.

Conditions: Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA