



# Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: [info@creativebiolabs.com](mailto:info@creativebiolabs.com)

**Product Name:** Malignant melanoma and normal adjacent tissue array  
**Catalog Number:** TMA- 120  
**Total Cases//Cores:** 50/100  
**Storage Conditions:** Cool Dry  
**Array Panel:** Duplicated cores per case, 45 cases of malignant melanoma and 5 cases normal tissue  
**Array Core Diameter:** 1.0 mm  
**Section Thickness:** 5 micrometers

## TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9	10
	A										
	B										
	C										
	D										
	E										
	F										
	G										
	H										
	I										
	J										

## PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	53	M	Skin	Malignant melanoma of right sole	-	
A2	53	M	Skin	Malignant melanoma of right sole	-	
A3	47	F	Skin	Malignant melanoma of right upper arm	-	
A4	47	F	Skin	Malignant melanoma of right upper arm	-	
A5	53	F	Skin	Malignant melanoma of right big toe	-	
A6	53	F	Skin	Malignant melanoma of right big toe	-	
A7	57	M	Skin	Malignant melanoma of right cheek	-	
A8	57	M	Skin	Malignant melanoma of right cheek	-	
A9	72	F	Skin	Malignant melanoma of vulva	-	
A10	72	F	Skin	Malignant melanoma of vulva	-	
B1	45	F	Skin	Malignant melanoma of vulva	-	
B2	45	F	Skin	Malignant melanoma of vulva	-	
B3	48	M	Skin	Malignant melanoma of left thumb	-	
B4	48	M	Skin	Malignant melanoma of left thumb	-	
B5	45	M	Skin	Malignant melanoma of perianal	-	
B6	45	M	Skin	Malignant melanoma of perianal	-	
B7	65	M	Skin	Malignant melanoma of scalp	-	
B8	65	M	Skin	Malignant melanoma of scalp	-	
B9	81	F	Skin	Malignant melanoma of left arm	-	
B10	81	F	Skin	Malignant melanoma of left arm	-	
C1	66	M	Skin	Malignant melanoma of right thigh	-	
C2	66	M	Skin	Malignant melanoma of right thigh	-	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828  
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



# Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: [info@creativebiolabs.com](mailto:info@creativebiolabs.com)

C3	35	M	Skin	Malignant melanoma of left dorsum of foot	-	
C4	35	M	Skin	Malignant melanoma of left dorsum of foot	-	
C5	51	M	Skin	Malignant melanoma of left oexter	-	
C6	51	M	Skin	Malignant melanoma of left oexter	-	
C7	71	M	Skin	Malignant melanoma of thigh	-	
C8	71	M	Skin	Malignant melanoma of thigh	-	
C9	36	M	Skin	Malignant melanoma of back	-	
C10	36	M	Skin	Malignant melanoma of back	-	
D1	49	M	Skin	Malignant melanoma of shoulder	-	
D2	49	M	Skin	Malignant melanoma of shoulder	-	
D3	52	M	Skin	Malignant melanoma of right heel	-	
D4	52	M	Skin	Malignant melanoma of right heel	-	
D5	52	F	Skin	Malignant melanoma of cunnus	-	
D6	52	F	Skin	Malignant melanoma of cunnus	-	
D7	38	F	Skin	Malignant melanoma of right toe	-	
D8	38	F	Skin	Malignant melanoma of right toe	-	
D9	25	M	Skin	Malignant melanoma of heel	-	
D10	25	M	Skin	Malignant melanoma of heel	-	
E1	66	F	Skin	Malignant melanoma of chest wall	-	
E2	66	F	Skin	Malignant melanoma of chest wall	-	
E3	62	F	Skin	Malignant melanoma of right thumb	-	
E4	62	F	Skin	Malignant melanoma of right thumb	-	
E5	77	F	Skin	Malignant melanoma of left sole	-	
E6	77	F	Skin	Malignant melanoma of left sole	-	
E7	40	M	Skin	Malignant melanoma of right chest wall	-	
E8	40	M	Skin	Malignant melanoma of right chest wall	-	
E9	62	M	Skin	Malignant melanoma of heel	-	
E10	62	M	Skin	Malignant melanoma of heel	-	
F1	47	F	Skin	Malignant melanoma of perianal	-	
F2	47	F	Skin	Malignant melanoma of perianal	-	
F3	38	M	Skin	Malignant melanoma of right abdominal wall	-	
F4	38	M	Skin	Malignant melanoma of right abdominal wall	-	
F5	58	F	Skin	Malignant melanoma of left buttock	-	
F6	58	F	Skin	Malignant melanoma of left buttock	-	
F7	42	F	Skin	Malignant melanoma of left thumb	-	
F8	42	F	Skin	Malignant melanoma of left thumb	-	
F9	54	F	Skin	Malignant melanoma of left sole	-	
F10	54	F	Skin	Malignant melanoma of left sole	-	
G1	57	M	Skin	Malignant melanoma of shoulder	-	
G2	57	M	Skin	Malignant melanoma of shoulder	-	
G3	88	F	Skin	Malignant melanoma of left sole	-	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828  
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



# Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: [info@creativebiolabs.com](mailto:info@creativebiolabs.com)

G4	88	F	Skin	Malignant melanoma of left sole	-	
G5	59	M	Skin	Malignant melanoma of left leg	-	
G6	59	M	Skin	Malignant melanoma of left leg	-	
G7	64	M	Skin	Malignant melanoma of left shoulder	-	
G8	64	M	Skin	Malignant melanoma of left shoulder	-	
G9	64	M	Skin	Malignant melanoma of sole	-	
G10	64	M	Skin	Malignant melanoma of sole	-	
H1	42	M	Skin	Malignant melanoma of left heel	-	
H2	42	M	Skin	Malignant melanoma of left heel	-	
H3	35	F	Skin	Malignant melanoma of waist	-	
H4	35	F	Skin	Malignant melanoma of waist	-	
H5	27	M	Skin	Little malignant melanoma of right thigh	-	
H6	27	M	Skin	Little malignant melanoma of right thigh	-	
H7	63	F	Skin	Malignant melanoma of right heel	-	
H8	63	F	Skin	Malignant melanoma of right heel	-	
H9	38	F	Skin	Little malignant melanoma of left upper arm	-	
H10	38	F	Skin	Little malignant melanoma of left upper arm	-	
I1	72	F	Skin	Malignant melanoma of vulva	-	
I2	72	F	Skin	Malignant melanoma of vulva	-	
I3	57	M	Skin	Malignant melanoma of right cheek	-	
I4	57	M	Skin	Malignant melanoma of right cheek	-	
I5	53	F	Skin	Malignant melanoma of right big toe	-	
I6	53	F	Skin	Malignant melanoma of right big toe	-	
I7	47	F	Skin	Malignant melanoma of right upper arm		
I8	47	F	Skin	Malignant melanoma of right upper arm	-	
I9	53	M	Skin	Malignant melanoma of right sole	-	
I10	53	M	Skin	Malignant melanoma of right sole	-	
J1	47	M	Skin	Skin tissue of wrist	-	
J2	47	M	Skin	Skin tissue of wrist	-	
J3	25	M	Skin	Skin tissue	-	
J4	25	M	Skin	Skin tissue	-	
J5	35	M	Skin	Skin tissue	-	
J6	35	M	Skin	Skin tissue	-	
J7	50	M	Skin	Skin tissue	-	
J8	50	M	Skin	Skin tissue	-	
J9	22	M	Skin	Skin tissue	-	
J10	22	M	Skin	Skin tissue	-	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828  
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



# Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: [info@creativebiolabs.com](mailto:info@creativebiolabs.com)

---

**For Research Use Only; Not For Diagnostic or Therapeutic Use.**

**Conditions:** Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

---

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828  
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA