


**Product Name:** Malignant melanoma with normal skin tissue array  
**Catalog Number:** TMA-1200  
**Total Cases//Cores:** 60/60  
**Storage Conditions:** Cool Dry  
**Array Panel:** Single core per case, 30 cases of carcinoma and 30 cases of adjacent tissue  
**Array Core Diameter:** 1.5 mm  
**Section Thickness:** 5 micrometers

## TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9	10
	A										
	B										
	C										
	D										
	E										
	F										

## PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	45	F	Skin	Malignant melanoma of vulva	-	
A2	72	F	Skin	Malignant melanoma of vulva	-	
A3	42	F	Skin	Malignant melanoma of left thumb	-	
A4	88	F	Skin	Malignant melanoma of left sole	-	
A5	58	F	Skin	Malignant melanoma of left buttock	-	
A6	54	F	Skin	Malignant melanoma of left sole	-	
A7	47	F	Skin	Malignant melanoma of crissum	-	
A8	65	M	Skin	Malignant melanoma of scalp	-	
A9	45	M	Skin	Malignant melanoma of crissum	-	
A10	52	M	Skin	Malignant melanoma of right heel	-	
B1	25	M	Skin	Malignant melanoma of heel	-	
B2	49	M	Skin	Malignant melanoma of shoulder	-	
B3	36	M	Skin	Malignant melanoma of back	-	
B4	71	M	Skin	Malignant melanoma of thigh	-	
B5	52	F	Skin	Malignant melanoma of vulva	-	
B6	66	F	Skin	Malignant melanoma of chest wall	-	
B7	77	F	Skin	Malignant melanoma of left heel	-	
B8	66	M	Skin	Malignant melanoma of right thigh	-	
B9	81	F	Skin	Malignant melanoma of left arm	-	
B10	38	F	Skin	Malignant melanoma of right toes	-	
C1	62	F	Skin	Malignant melanoma of right thumb	-	
C2	38	M	Skin	A little malignant melanoma of right abdominal wall	-	
C3	57	M	Skin	Malignant melanoma of right cheek	-	
C4	47	F	Skin	Malignant melanoma of right upper arm	-	
C5	57	M	Skin	Malignant melanoma of left scapula	-	



# Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: [info@creativebiolabs.com](mailto:info@creativebiolabs.com)

C6	16	F	Skin	Malignant melanoma of back	-	
C7	51	M	Skin	Malignant melanoma of left armpit	-	
C8	63	F	Skin	Malignant melanoma of right heel	-	
C9	74	M	Skin	Malignant melanoma of left foot	-	
C10	27	M	Skin	Malignant melanoma of right thigh	-	
D1	49	F	Skin	Skin tissue of right breast	-	
D2	46	F	Skin	Skin tissue of left breast	-	
D3	50	F	Skin	Skin tissue of right breast	-	
D4	49	F	Skin	Skin tissue of right breast	-	
D5	45	F	Skin	Skin tissue of left breast	-	
D6	42	F	Skin	Skin tissue of left breast	-	
D7	16	M	Skin	Skin tissue of right shoulder	-	
D8	40	F	Skin	Skin tissue of right elbow	-	
D9	43	F	Skin	Skin tissue of left breast	-	
D10	69	F	Skin	Skin tissue of left breast	-	
E1	36	M	Skin	Skin tissue of back	-	
E2	49	F	Skin	Skin tissue of breast	-	
E3	42	F	Skin	Skin tissue of left breast	-	
E4	58	F	Skin	Skin tissue of left breast	-	
E5	70	F	Skin	Skin tissue of left breast	-	
E6	75	F	Skin	Skin tissue of right hand	-	
E7	40	F	Skin	Skin tissue of breast	-	
E8	46	F	Skin	Skin tissue of left breast	-	
E9	73	M	Skin	A little epidermis and keratinized material of left-heel	-	
E10	44	F	Skin	Skin tissue of breast	-	
F1	51	M	Skin	Skin tissue of left upper abdominal wall	-	
F2	45	M	Skin	Chronic inflammation of skin tissue of scalp	-	
F3	58	F	Skin	Skin tissue of right breast	-	
F4	49	F	Skin	Skin tissue	-	
F5	49	F	Skin	Skin tissue	-	
F6	30	M	Skin	Skin tissue of foot	-	
F7	35	M	Skin	Skin tissue	-	
F8	22	M	Skin	Skin tissue	-	
F9	47	M	Skin	Skin tissue of wrist	-	
F10	35	M	Skin	Skin tissue of foot	-	

## For Research Use Only; Not For Diagnostic or Therapeutic Use.

**Conditions:** Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828  
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA