




Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

Product Name: Colon adenocarcinoma tissue array
Catalog Number: TMA-039
Total Cases/Cores: 23/63
Storage Conditions: Cool Dry
Array Panel: Triplicated cores per case, 20 cases of carcinoma and 3 cases of adjacent tissue
Array Core Diameter: 1.5 mm
Section Thickness: 5 micrometers

TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9
	A									
	B									
	C									
	D									
	E									
	F									
	G									

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	62	F	Colon	Adenocarcinoma	Well	
A2	62	F	Colon	Adenocarcinoma	Well	
A3	62	F	Colon	Adenocarcinoma	Well	
A4	26	M	Colon	Adenocarcinoma	Well	
A5	26	M	Colon	Adenocarcinoma	Well	
A6	26	M	Colon	Adenocarcinoma	Well	
A7	51	F	Colon	Mucous adenocarcinoma	Well	
A8	51	F	Colon	Adenocarcinoma	Well	
A9	51	F	Colon	Mucous adenocarcinoma	Well	
B1	55	M	Colon	Adenocarcinoma	Well	
B2	55	M	Colon	Adenocarcinoma	Well	
B3	55	M	Colon	Adenocarcinoma	Well	
B4	59	M	Colon	Adenocarcinoma	Moderately	
B5	59	M	Colon	Adenocarcinoma	Well	
B6	59	M	Colon	Adenocarcinoma	Well	
B7	57	F	Colon	Adenocarcinoma	Well	
B8	57	F	Colon	Adenocarcinoma	Well	
B9	57	F	Colon	Adenocarcinoma	Well	
C1	86	M	Colon	Adenocarcinoma	Moderately	
C2	86	M	Colon	Adenocarcinoma	Moderately	
C3	86	M	Colon	Adenocarcinoma	Moderately	
C4	53	F	Colon	Adenocarcinoma	Well	
C5	53	F	Colon	Adenocarcinoma	Well	
C6	53	F	Colon	Adenocarcinoma	Well	
C7	20	M	Colon	Adenocarcinoma	Poorly	
C8	20	M	Colon	Adenocarcinoma	Poorly	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C9	20	M	Colon	Adenocarcinoma	Poorly	
D1	75	M	Colon	Adenocarcinoma	Well	
D2	75	M	Colon	Adenocarcinoma	Well	
D3	75	M	Colon	Adenocarcinoma	Well	
D4	81	M	Colon	Adenocarcinoma	Well	
D5	81	M	Colon	Adenocarcinoma	Well	
D6	81	M	Colon	Adenocarcinoma	Well	
D7	72	M	Colon	Adenocarcinoma	Well	
D8	72	M	Colon	Adenocarcinoma	Well	
D9	72	M	Colon	Adenocarcinoma	Well	
E1	48	M	Colon	Adenocarcinoma	Poorly	
E2	48	M	Colon	Adenocarcinoma	Poorly	
E3	48	M	Colon	Adenocarcinoma	Poorly	
E4	59	F	Colon	Adenocarcinoma	Well	
E5	59	F	Colon	Adenocarcinoma	Well	
E6	59	F	Colon	Adenocarcinoma	Well	
E7	44	M	Colon	Adenocarcinoma	Well	
E8	44	M	Colon	Adenocarcinoma	Well	
E9	44	M	Colon	Adenocarcinoma	Well	
F1	67	F	Colon	Adenocarcinoma	Well	
F2	67	F	Colon	Adenocarcinoma	Well	
F3	67	F	Colon	Adenocarcinoma	Well	
F4	64	F	Colon	Adenocarcinoma	Well	
F5	64	F	Colon	Adenocarcinoma	Well	
F6	64	F	Colon	Adenocarcinoma	Well	
F7	60	M	Colon	Adenocarcinoma	Well	
F8	60	M	Colon	Adenocarcinoma	Well	
F9	60	M	Colon	Adenocarcinoma	Well	
G1	47	M	Colon	Adenocarcinoma	Well	
G2	47	M	Colon	Adenocarcinoma	Well	
G3	47	M	Colon	Adenocarcinoma	Well	
G4	62	M	Colon	Adenocarcinoma	Well	
G5	62	M	Colon	Adenocarcinoma	Well	
G6	62	M	Colon	Adenocarcinoma	Well	
G7	69	M	Small intestine	Mucosa of small intestine	-	
G8	56	F	Colon	Mucosa	-	
G9	70	M	Colon	Mucosa	-	

For Research Use Only; Not For Diagnostic or Therapeutic Use.

Conditions: Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA