


Product Name: Multiple ovarian carcinoma tissue array
Catalog Number: TMA-461
Total Cases//Cores: 70/70
Storage Conditions: Cool Dry
Array Panel: Single core per case, 70 cases of ovarian carcinoma
Array Core Diameter: 1.5 mm
Section Thickness: 5 micrometers

TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9	10
	A										
	B										
	C										
	D										
	E										
	F										
	G										

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	40	F	Ovary	Clear cell carcinoma	-	
A2	57	F	Ovary	Serous papillary carcinoma	Moderate	
A3	48	F	Ovary	Clear cell carcinoma	-	
A4	57	F	Ovary	Serous papillary carcinoma	Moderate	
A5	43	F	Ovary	Serous papillary carcinoma	Poorly	
A6	54	F	Ovary	Serous papillary carcinoma	Well	
A7	63	F	Ovary	Serous papillary carcinoma	Poorly	
A8	46	F	Ovary	Serous papillary carcinoma	Poorly	
A9	54	F	Ovary	Serous papillary carcinoma	Poorly	
A10	56	F	Ovary	Hyperplastic fibrous tissue	-	
B1	44	F	Ovary	Granular cell tumor	-	
B2	49	F	Ovary	Serous papillary carcinoma	Poorly	
B3	18	F	Ovary	Immature teratoma	-	
B4	15	F	Ovary	Endodermal sinus carcinoma	-	
B5	38	F	Ovary	Metastatic adenocarcinoma	Poorly	
B6	39	F	Ovary	Serous papillary carcinoma	Poorly	
B7	24	F	Ovary	Endodermal sinus carcinoma	-	
B8	42	F	Ovary	Serous papillary carcinoma	Poorly	
B9	50	F	Ovary	Serous papillary carcinoma	Poorly	
B10	49	F	Ovary	Serous papillary carcinoma	Poorly	
C1	62	F	Ovary	Serous papillary carcinoma	Poorly	
C2	53	F	Ovary	Mucinous papillary carcinoma	Well	
C3	38	F	Ovary	Metastatic adenocarcinoma	Poorly	
C4	43	F	Ovary	Clear cell carcinoma	-	
C5	26	F	Ovary	Serous papillary carcinoma	Well	
C6	47	F	Ovary	Serous papillary carcinoma	Well	
C7	62	F	Ovary	Squamous cell carcinoma	Poorly	



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C8	35	F	Ovary	Dysgerminoma	-	
C9	41	F	Ovary	Dysgerminoma	-	
C10	47	F	Ovary	Serous papillary carcinoma	Poorly	
D1	42	F	Ovary	Clear cell carcinoma	-	
D2	39	F	Ovary	Metastatic adenocarcinoma	Well	
D3	66	F	Ovary	Metastatic adenocarcinoma	Poorly	
D4	48	F	Ovary	Malignant follicular theca cytoma	-	
D5	51	F	Ovary	Serous papillary carcinoma	Well	
D6	33	F	Ovary	Metastatic signet-ring cell carcinoma	-	
D7	18	F	Ovary	Mixed germ cell tumor	-	
D8	40	F	Ovary	Metastatic signet-ring cell carcinoma	-	
D9	43	F	Ovary	Granular cell tumor	-	
D10	55	F	Ovary	Serous papillary carcinoma	Poorly	
E1	46	F	Ovary	Serous papillary carcinoma	Poorly	
E2	57	F	Ovary	Serous papillary carcinoma	Poorly	
E3	75	F	Ovary	Serous papillary carcinoma	Poorly	
E4	69	F	Ovary	Little serous papillary carcinoma	Poorly	
E5	30	F	Ovary	Serous papillary carcinoma	Poorly	
E6	42	F	Ovary	Serous papillary carcinoma	Moderately	
E7	48	F	Ovary	Clear cell carcinoma	-	
E8	22	F	Ovary	Serous papillary carcinoma	Well	
E9	50	F	Ovary	Clear cell carcinoma	-	
E10	32	F	Ovary	Serous papillary carcinoma	Moderately	
F1	48	F	Ovary	Serous papillary carcinoma	Moderately	
F2	50	F	Ovary	Serous papillary carcinoma	Moderately	
F3	65	F	Ovary	Serous papillary carcinoma	Poorly	
F4	38	F	Ovary	Serous papillary carcinoma	Poorly	
F5	31	F	Ovary	Metastatic adenocarcinoma	Poorly	
F6	55	F	Ovary	Metastatic adenocarcinoma	Poorly	
F7	51	F	Ovary	Serous papillary carcinoma	Poorly	
F8	65	F	Ovary	Serous papillary carcinoma	Moderately	
F9	26	F	Ovary	Serous papillary carcinoma	Moderately	
F10	55	F	Ovary	Serous papillary carcinoma	Moderately	
G1	49	F	Ovary	Serous papillary carcinoma	Poorly	
G2	48	F	Ovary	Metastatic signet-ring cell carcinoma	-	
G3	46	F	Ovary	Serous papillary carcinoma	Poorly	
G4	63	F	Ovary	Serous papillary carcinoma	Poorly	
G5	37	F	Ovary	Serous papillary carcinoma	Moderately	
G6	35	F	Ovary	Little malignant tumor	-	
G7	12	F	Ovary	Dysgerminoma	-	
G8	55	F	Ovary	Serous papillary carcinoma	Moderately	
G9	20	F	Ovary	Little malignant tumor cell	-	
G10	55	F	Ovary	Serous papillary carcinoma	Moderately	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA