


Product Name: Cervix carcinoma tissue array
Catalog Number: TMA-501
Total Cases//Cores: 54/54
Storage Conditions: Cool Dry
Array Panel: Single core per case, 40 cases of carcinoma, 10 cases of adjacent tissue and 4 cases normal tissue
Array Core Diameter: 1.5 mm
Section Thickness: 5 micrometers

TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9	10
	A										
	B										
	C										
	D										
	E										
	F										

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	40	F	Uterine cervix	Squamous cell carcinoma	Well	
A2	56	F	Uterine cervix	Squamous cell carcinoma	Well	
A3	47	F	Uterine cervix	Squamous cell carcinoma	Well	
A4	65	F	Uterine cervix	Squamous cell carcinoma	Well	
A5	53	F	Uterine cervix	Squamous cell carcinoma	Well	
A6	69	F	Uterine cervix	Squamous cell carcinoma	Well	
A7	38	F	Uterine cervix	Fibrous tissue	-	
A8	48	F	Uterine cervix	A little squamous cell carcinoma	Well	
A9	48	F	Uterine cervix	Squamous cell carcinoma	Well	
A10	43	F	Uterine cervix	Squamous cell carcinoma	Well	
B1	52	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B2	36	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B3	63	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B4	27	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B5	47	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B6	68	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B7	53	F	Uterine cervix	Squamous cell carcinoma	Well	
B8	48	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B9	76	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B10	37	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C1	27	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C2	38	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C3	36	F	Uterine cervix	Squamous cell carcinoma	Poorly	



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C4	39	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C5	37	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C6	27	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C7	58	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C8	27	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C9	51	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C10	32	F	Uterine cervix	Squamous cell carcinoma	Poorly	
D1	48	F	Uterine cervix	Adenosquamous carcinoma	-	
D2	43	F	Uterine cervix	Adenosquamous carcinoma	-	
D3	38	F	Uterine cervix	Adenosquamous carcinoma	-	
D4	54	F	Uterine cervix	Adenosquamous carcinoma	-	
D5	45	F	Uterine cervix	Adenosquamous carcinoma	-	
D6	42	F	Uterine cervix	Adenosquamous carcinoma	-	
D7	43	F	Uterine cervix	Adenosquamous carcinoma	-	
D8	54	F	Uterine cervix	Adenosquamous carcinoma	-	
D9	61	F	Uterine cervix	Adenosquamous carcinoma	-	
D10	42	F	Uterine cervix	Adenosquamous carcinoma	-	
E1	40	F	Uterine cervix	Cervical tissue	-	
E2	41	F	Uterine cervix	Chronic inflammation with CIN II	-	
E3	46	F	Uterine cervix	Cervical tissue	-	
E4	32	F	Uterine cervix	Little chronic inflammation of cervix	-	
E5	41	F	Uterine cervix	Squamous cell carcinoma	Well	
E6	49	F	Uterine cervix	Cervical canals tissue	-	
E7	37	F	Uterine cervix	Cervical tissue	-	
E8	35	F	Uterine cervix	Cervical canals tissue	-	
E9	32	F	Uterine cervix	Chronic inflammation of cervix	-	
E10	40	F	Uterine cervix	Cervical canals tissue	-	
F1	38	F	Uterine cervix	Little cervical tissue	-	
F2	21	F	Uterine cervix	Fibrous tissue	-	
F3	21	F	Uterine cervix	Cervical tissue	-	
F4	15	F	Uterine cervix	Little cervical tissue	-	

For Research Use Only; Not For Diagnostic or Therapeutic Use.

Conditions: Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA