


**Product Name:** Cervix carcinoma and adjacent tissue array  
**Catalog Number:** TMA-502  
**Total Cases//Cores:** 50/50  
**Storage Conditions:** Cool Dry  
**Array Panel:** Single core per case, 46 cases of carcinoma, 4 cases of adjacent tissue  
**Array Core Diameter:** 1.5 mm  
**Section Thickness:** 5 micrometers

### TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9	10
	A										
	B										
	C										
	D										
	E										

### PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	56	F	Uterine cervix	Squamous cell carcinoma	Well	
A2	49	F	Uterine cervix	Squamous cell carcinoma	Well	
A3	60	F	Uterine cervix	Squamous cell carcinoma	Well	
A4	54	F	Uterine cervix	Squamous cell carcinoma	Well	
A5	32	F	Uterine cervix	Squamous cell carcinoma	Well	
A6	40	F	Uterine cervix	Squamous cell carcinoma	Moderately	
A7	65	F	Uterine cervix	Squamous cell carcinoma	Well	
A8	53	F	Uterine cervix	Squamous cell carcinoma	Well	
A9	48	F	Uterine cervix	Squamous cell carcinoma	Well	
A10	52	F	Uterine cervix	Squamous cell carcinoma	Well	
B1	36	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B2	37	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B3	53	F	Uterine cervix	Squamous cell carcinoma	Well	
B4	36	F	Uterine cervix	A little squamous cell carcinoma	Poorly	
B5	48	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B6	46	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B7	39	F	Uterine cervix	Squamous cell carcinoma	Moderately	
B8	36	F	Uterine cervix	Squamous cell carcinoma	Well	
B9	25	F	Uterine cervix	Fibrous and blood vessel tissue	-	
B10	69	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C1	58	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C2	27	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C3	68	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C4	37	F	Uterine cervix	Squamous cell carcinoma	Poorly	



# Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: [info@creativebiolabs.com](mailto:info@creativebiolabs.com)

C5	63	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C6	43	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C7	47	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C8	38	F	Uterine cervix	Squamous cell carcinoma	Moderately	
C9	27	F	Uterine cervix	Squamous cell carcinoma	Poorly	
C10	49	F	Uterine cervix	Squamous cell carcinoma	Poorly	
D1	62	F	Uterine cervix	Squamous cell carcinoma	Moderately	
D2	36	F	Uterine cervix	Squamous cell carcinoma	Moderately	
D3	76	F	Uterine cervix	Squamous cell carcinoma	Poorly	
D4	37	F	Uterine cervix	Squamous cell carcinoma	Moderately	
D5	27	F	Uterine cervix	Squamous cell carcinoma	Poorly	
D6	40	F	Uterine cervix	Squamous cell carcinoma	Moderately	
D7	27	F	Uterine cervix	Squamous cell carcinoma	Poorly	
D8	39	F	Uterine cervix	Squamous cell carcinoma	Moderately	
D9	39	F	Uterine cervix	Squamous cell carcinoma	Poorly	
D10	48	F	Uterine cervix	Squamous cell carcinoma	Poorly	
F1	42	F	Uterine cervix	Mucous adenocarcinoma	Well	
F2	44	F	Uterine cervix	Adenocarcinoma	Moderately	
F3	32	F	Uterine cervix	Fibrous and blood vessel tissue	-	
F4	43	F	Uterine cervix	Adenosquamous carcinoma	-	
F5	64	F	Uterine cervix	Adenosquamous carcinoma	-	
F6	38	F	Uterine cervix	Adenosquamous carcinoma	-	
F7	21	F	Uterine cervix	A little cervix tissue	-	
F8	18	F	Uterine cervix	Cervix tissue	-	
F9	15	F	Uterine cervix	Cervix tissue	-	
F10	21	F	Uterine cervix	Cervical canals tissue	-	

## For Research Use Only; Not For Diagnostic or Therapeutic Use.

**Conditions:** Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828  
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA