




Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

Product Name: Endometrial carcinoma (multi-tissue combined panel) tissue array
Catalog Number: TMA-511
Total Cases/Cores: 63/63
Storage Conditions: Cool Dry
Array Panel: Single core per case, 60 cases of carcinoma and 3 cases of adjacent tissue
Array Core Diameter: 1.5 mm
Section Thickness: 5 micrometers

TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9
	A									
	B									
	C									
	D									
	E									
	F									
	G									

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	61	F	Uterus	A little carcinoma tissue	-	
A2	52	F	Uterus	Endometrial adenocarcinoma	Moderately	
A3	64	F	Uterus	Endometrial adenocarcinoma	Well	
A4	52	F	Uterus	Endometrial adenocarcinoma	Well	
A5	51	F	Uterus	Endometrial adenocarcinoma	Well	
A6	58	F	Uterus	Endometrial adenocarcinoma	Poorly	
A7	44	F	Uterus	Endometrial adenocarcinoma	Moderately	
A8	54	F	Uterus	Endometrial adenocarcinoma	Moderately	
A9	38	F	Uterus	Endometrial adenocarcinoma	Well	
B1	48	F	Uterus	Endometrial adenocarcinoma	Well	
B2	55	F	Uterus	Invasive hydatidiform mole	-	
B3	60	F	Uterus	Endometrial adenocarcinoma	Moderately	
B4	63	F	Uterus	Endometrial adenocarcinoma	Moderately	
B5	53	F	Uterus	Endometrial adenocarcinoma	Moderately	
B6	53	F	Uterus	Endometrial adenocarcinoma	Well	
B7	66	F	Uterus	Endometrial adenocarcinoma	Moderately	
B8	56	F	Uterus	Endometrial adenocarcinoma	Moderately	
B9	49	F	Uterus	Endometrial adenocarcinoma	Moderately	
C1	54	F	Uterus	Endometrial adenocarcinoma	Moderately	
C2	58	F	Uterus	Endometrial adenocarcinoma	Well	
C3	62	F	Uterus	Endometrial adenocarcinoma	Well	
C4	57	F	Uterus	Endometrial adenocarcinoma	Moderately	
C5	51	F	Uterus	Endometrial adenocarcinoma	Moderately	
C6	45	F	Uterus	Endometrial adenocarcinoma	Moderately	
C7	55	F	Uterus	Endometrial adenocarcinoma	Well	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C8	59	F	Uterus	Endometrial adenocarcinoma	Moderately	
C9	51	F	Uterus	Endometrial adenocarcinoma	Moderately	
D1	55	F	Uterus	Endometrial adenocarcinoma	Well	
D2	50	F	Uterus	Endometrial adenocarcinoma	Well	
D3	80	F	Uterus	Endometrial adenocarcinoma	Well	
D4	65	F	Uterus	Adenosquamous carcinoma	-	
D5	62	F	Uterus	Endometrial adenocarcinoma	Moderately	
D6	58	F	Uterus	Endometrial adenocarcinoma	Moderately	
D7	57	F	Uterus	Endometrial adenocarcinoma	Well	
D8	51	F	Uterus	Endometrial adenocarcinoma	Moderately	
D9	57	F	Uterus	Endometrial adenocarcinoma	Well	
E1	67	F	Uterus	Endometrial adenocarcinoma	Moderately	
E2	45	F	Uterus	Endometrial adenocarcinoma	Moderately	
E3	62	F	Uterus	Endometrial adenocarcinoma	Poorly	
E4	52	F	Uterus	Endometrial adenocarcinoma	Moderately	
E5	55	F	Uterus	Endometrial adenocarcinoma	Moderately	
E6	72	F	Uterus	Endometrial adenocarcinoma	Moderately	
E7	62	F	Uterus	Endometrial adenocarcinoma	Moderately	
E8	49	F	Uterus	Adenosquamous carcinoma	-	
E9	51	F	Uterus	Endometrial adenocarcinoma	Moderately	
F1	62	F	Uterus	Endometrial adenocarcinoma	Well	
F2	53	F	Uterus	Endometrial adenocarcinoma	Poorly	
F3	48	F	Uterus	Endometrial adenocarcinoma	Poorly	
F4	39	F	Uterus	Endometrial adenocarcinoma	Moderately	
F5	60	F	Uterus	Endometrial papillary adenocarcinoma	Well	
F6	68	F	Uterus	Endometrial adenocarcinoma	Well	
F7	47	F	Uterus	Endometrial adenocarcinoma	Moderately	
F8	52	F	Uterus	Endometrial adenocarcinoma	Moderately	
F9	56	F	Uterus	Endometrial adenocarcinoma	Moderately	
G1	63	F	Uterus	Endometrial adenocarcinoma	Poorly	
G2	54	F	Uterus	Endometrial adenocarcinoma	Moderately	
G3	72	F	Uterus	Endometrial adenocarcinoma	Moderately	
G4	50	F	Uterus	Endometrial adenocarcinoma	Moderately	
G5	78	F	Uterus	Endometrial adenocarcinoma	Poorly	
G6	56	F	Uterus	Endometrial adenocarcinoma	Well	
G7	48	F	Uterus	Secretory endometrium	-	
G8	56	F	Uterus	Proliferative endometrium	-	
G9	61	F	Uterus	Smooth muscle, blood vessel	-	

For Research Use Only; Not For Diagnostic or Therapeutic Use.

Conditions: Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA