


Product Name: Breast carcinoma tissue array
Catalog Number: CBL-TMA-555
Total Cases/Cores: 63/63
Storage Conditions: Cool Dry
Array Panel: Single core per case, 63 cases of carcinoma
Array Core Diameter: 1.5 mm
Section Thickness: 5 micrometers

TISSUE CORE PANEL DISPLAY ON SLIDE

		1	2	3	4	5	6	7	8	9
	A									
	B									
	C									
	D									
	E									
	F									
	G									

PRODUCT SPECIFICATION SHEET

No.	Age	Sex	Organ	Pathology diagnosis	Differentiation	Note
A1	37	F	Breast	Breast tissue	-	
A2	40	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
A3	40	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
A4	48	F	Breast	Atypical medullary carcinoma	-	
A5	31	F	Breast	Atypical medullary carcinoma	-	
A6	72	F	Breast	Atypical medullary carcinoma	-	
A7	47	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
A8	57	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
A9	54	F	Breast	Infiltrating lobular carcinoma	-	
B1	50	M	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
B2	40	F	Breast	Atypical medullary carcinoma	-	
B3	50	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
B4	41	F	Breast	Nonspecific infiltrating ductal carcinoma	Well-Moderately	
B5	32	F	Breast	Atypical medullary carcinoma	-	
B6	43	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately-Poorly	
B7	70	F	Breast	A little nonspecific infiltrating ductal carcinoma	-	
B8	47	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
B9	45	F	Breast	Complex lobule and ductal carcinoma	-	
C1	49	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
C2	52	F	Breast	Infiltrating lobular carcinoma	-	
C3	48	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
C4	38	F	Breast	Infiltrating lobular carcinoma	-	
C5	68	F	Breast	Nonspecific infiltrating ductal carcinoma	Well-Moderately	



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

C6	60	F	Breast	Infiltrating lobular carcinoma	-	
C7	35	F	Breast	Infiltrating lobular carcinoma	-	
C8	53	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
C9	38	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D1	55	M	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D2	54	F	Breast	A little atypical medullary carcinoma	-	
D3	46	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D4	48	F	Breast	Nonspecific infiltrating ductal carcinoma	Well- Moderately	
D5	49	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
D6	51	F	Breast	Complex lobule and ductal carcinoma	-	
D7	32	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
D8	52	F	Breast	A little nonspecific infiltrating ductal carcinoma	-	
D9	47	F	Breast	Nonspecific infiltrating ductal carcinoma	Well- Moderately	
E1	50	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
E2	55	F	Breast	Atypical medullary carcinoma	-	
E3	42	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E4	53	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E5	70	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
E6	28	F	Breast	Infiltrating lobular carcinoma	-	
E7	44	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E8	52	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
E9	45	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F1	62	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
F2	48	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F3	38	F	Breast	Infiltrating lobular carcinoma	-	
F4	38	F	Breast	Complex lobule and ductal carcinoma	-	
F5	42	F	Breast	Nonspecific infiltrating ductal carcinoma	Well	
F6	48	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F7	52	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
F8	45	F	Breast	Nonspecific infiltrating ductal carcinoma	Poorly	
F9	49	F	Breast	Atypical medullary carcinoma	-	
G1	65	F	Breast	Nonspecific infiltrating ductal carcinoma	Well- Moderately	
G2	55	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G3	39	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G4	60	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G5	59	F	Breast	Nonspecific infiltrating ductal carcinoma	Well- Moderately	
G6	61	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
G7	47	F	Breast	Infiltrating lobular carcinoma	-	
G8	49	F	Breast	Complex lobule and ductal carcinoma	-	

<http://www.creativebiolabs.com>

Tel: (631)-871-5806, Fax: (631)-614-7828
21 Brookhaven Blvd. Port Jefferson Station, NY 11776, USA



Product Information Sheet

Creative Dynamics Inc Phone: (631) 871-5806 Fax: (631) 614-7828 Email: info@creativebiolabs.com

G9	43	F	Breast	Nonspecific infiltrating ductal carcinoma	Moderately	
----	----	---	--------	---	------------	--

For Research Use Only; Not For Diagnostic or Therapeutic Use.

Conditions: Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of CD Inc is strictly prohibited.