

Product Information

YPD Plus Medium

Cat. No.: **HN-1124-TT14**

This product is for research use only and is not intended for diagnostic use.

Product Overview

YPD Plus contains extremely rich nutrients. Due to the special high-pressure sterilization, it can be sealed and temporarily stored at 4 °C or room temperature. To avoid pollution, please store at -20 ° C for long-term storage. When using YPD Plus Medium recovery and other methods to improve the transformation efficiency of yeast, PEG/LiAc transformation solution, Carrier DNA, DMSO, logarithmic yeast, and optimized plasmid DNA quality and concentration are usually added. Recovering with YPD Plus for 30-60 minutes can effectively improve transformation efficiency by over 60%. This product is suitable for *Saccharomyces cerevisiae* such as Y2HGold and Y1HGold, especially for library transformation that requires higher transformation efficiency.

Applications

Library construction

Size

50ml

Storage

4°C/-20°C

Recommended dosage

1.3g/L