

Hirudin

Recombinant (rHirudin)

Expressed in *Pichia. Pastoris*

Cat. No. CRP08123

Lot. No. (See product label)

PRODUCT INFORMATION

Description: Hirudin is a potent thrombin inhibitor originally derived from the medicinal leech. Unlike heparin, hirudin acts directly on thrombin, rather than through other clotting factors. They have a high binding affinity and specificity for thrombin. The mechanism of hirudin-thrombin binding appears to be unique. Recombinant hirudin variant is derived from yeast cell, the polypeptide containing 65 amino acid residues has a molecular weight of 6979.5 Da, which is identical to natural hirudin except for substitution of leucine for isoleucine at the N-terminal end of the molecule and the absence of a sulfate group on the tyrosine at position 63.

M. W. : 6979.5 Da

Amino-Acid Sequence: 65aa. non-glycosylated

CAS Number: 8001-27-2

Recombinant: Expressed in *Pichia. Pastoris*

Purity: >96% by SDS-PAGE and HPLC analyses.

Specific Activity: rHirudin is fully biologically active when compared to standard. Its specific activity is $\geq 1 \times 10^4$ ATU/mg.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Production: Yeast expression of a DNA sequence encoding the mature hirudin variant (Leu1,Thr2-63-desulfo hirudin).

Endotoxin: Less than 10EU/mg of rHirudin as determined by LAL method.

Formulation: Lyophilized from a 0.2 μ m filtered solution of 20mM PBS, pH 7.0, containing 2% mannitol.

Reconstitution: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at <-20°C. Further dilutions should be made in appropriate buffered solutions.

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C. Avoid repeated freeze/thaw cycles.

REFERENCES

1. Folkers PJM, Clore GM. et al. (1989). Solution structure of recombinant hirudin and the Lys-47-Glu mutant: a nuclear magnetic resonance and hybrid distance geometry-dynamical simulated annealing study. *Biochemistry* 28 (6): 2601-2617.
2. Haruyama H. and Wuthrich K. (1989). Conformation of recombinant desulfatohirudin in aqueous solution determined by nuclear magnetic resonance. *Biochemistry* 28 (10): 4301-4312.
3. Rydell TJ, Tulinsky A. et al. (1991). Refined structure of the Hirudin-Thrombin complex. *J. Mol. Biol.* 221 (2): 583-601.

FOR RESEARCH USE ONLY

2005-2008 Creative Biolabs. All rights reserved.