

Product Information

Saccharomyces cerevisiae Competent cell Preparation Kit

Cat. No.: **HN-0125-TT1**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Saccharomyces cerevisiae competent cell Preparation Kit is mainly used for Saccharomyces cerevisiae plasmid transformation. This kit includes: TE/LiAC yeast competent preparation solution (frozen type), PEG/LiAC yeast transformation solution, as well as Carrier DNA for transformation and a complimentary bottle of YPD Plus. The frozen state preparation solution and transformation solution in this reagent kit are developed through special processes, with well proportioned working solutions that can be directly used and have stable quality for long-term storage. The sensitive state prepared by this reagent kit can be stored at -80 °C for easy transformation and use in the later stage.

Applications

Yeast hybrid

Content

TE/LiAC yeast competent preparation solution (frozen type)

PEG/LiAC yeast transformation solution

Carrier DNA

YPD Plus

Size

100T/200T

Storage

Please refer to the manual.