

# Product Information

## One-step Yeast Transformation Kit

Cat. No.: **HN-0125-TT18**

This product is for research use only and is not intended for diagnostic use.

### Product Overview

This one-step yeast fast transformation kit has been developed for *Saccharomyces cerevisiae* through a special process, which does not require the preparation of competent cells. Only one step of heat shock is needed to complete the entire transformation process within 45 minutes, greatly reducing the transformation time. The process is simple and convenient, and suitable for the transformation of non library plasmids, as well as for the transformation of bait plaques such as Y1H [pAbAi Pro]. Its transformation rate can reach 60-85% of the classical yeast transformation method (Cat. A04LN-0125-TT4). The freshness of yeast spots is directly proportional to their transformation efficiency, and it is recommended to use fresh yeast spots as much as possible when selecting spots.

### Applications

Yeast transformation

### Content

YS1 transformation solution

YPD Plus

### Size

50T

### Storage

Please refer to the manual.