

Product Information

X- α -Gal

Cat. No.: **HN-0125-TT8**

This product is for research use only and is not intended for diagnostic use.

Product Overview

X- α -Gal is a chromogenic substrate of yeast galactosidase (MEL1) used for direct detection of yeast two hybrid activity in GAL4 system on culture medium. Through blue/white color screening, positive blue clones can be quickly and easily identified without the time-consuming β - galactosidase reporter gene testing.

Applications

Yeast hybrid

Content

X- α -Gal powder

Size

250mg/5×250mg/500mg/1g

Storage

-20°C

Protected from light