

Product Information

50%PEG Solution

Cat. No.: **HN-0125-TT9**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Polyethylene glycol (PEG) has many varieties, including PEG400, PEG1500, PEG3350, PEG4000, PEG6000, PEG8000, etc. The PEG series products are non-toxic, non irritating, slightly bitter in taste, have good water solubility, and have good compatibility with many organic components. They have excellent lubricity, moisturizing, dispersibility, adhesives, anti-static agents, and softeners.

This product is a sterile 50% PEG solution, commonly used to promote cell fusion or protoplast fusion and aid in the uptake of DNA by organisms (such as yeast) during transformation.

Applications

Yeast transformation

Content

50%PEG solution

Size

50ml

Storage

RT