

Product Information

SG/-His/-Leu/-Trp/-Ura Broth

Cat. No.: **HN-1124-TT156**

This product is for research use only and is not intended for diagnostic use.

Product Overview

SG type culture medium replaces glucose with galactose in SD type culture medium, with the same components as SD type culture medium except for the carbon source. For example, SG/-Ura Broth is a mixture of Minimum SD Base Gal and DO Supplement Ura.

We also have medium with Agar added, such as SG/-Ura With Agar. Compared with conventional SD/-Ura Broth, different carbon source culture media include:

SG/-Ura Broth (containing galactose), SF/-Ura Broth (containing fructose), SM/-Ura Broth (containing maltose), SS/-Ura Broth (containing sucrose), etc.

The culture medium with different carbon sources also includes different defect types, which can be selected according to experimental needs.

Applications

Yeast expression and hybrid

Size

10×0.5L

Storage

RT

Recommended dosage

28g/L