
Product Information

Yeast Specific Agar

Cat. No.: **HN-1124-TT4**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Agar powder is composed of agarose and agarose pectin, which is insoluble in cold water and easily soluble in boiling water, and slowly soluble in hot water. This is a high-purity agar powder. After the agar powder solution with a gel strength of 1.5% solidifies at 20 degrees Celsius, the maximum pressure that the surface can withstand without breaking is greater than 1300g/cm².

Applications

Yeast cultivation

Size

500g

Storage

RT

Recommended dosage

1.5-2% or 6-25g/L