

Product Information

Custom Yeast Hybrid Media Kit

Cat. No.: **HN-1224-TT20**

This product is for research use only and is not intended for diagnostic use.

Product Overview

During yeast hybrid experiments, different types of screening media are required based on the experimental objectives and protocols. We have selected 14 commonly used culture media for you to freely mix and match to meet different needs. The set includes three series of products: Rich Medium, SD Broth, and SD with Agar, which meet the needs of all culture media in your Y1H, Y2H, and Y3H processes and avoid unnecessary waste.

Applications

Yeast hybrid

Content

Custom component

Size

1Set

Storage

Please refer to the manual.