

Product Information

Y1HGold-pAbAi Yeast 1 Hybrid Library Screening Kit

Cat. No.: **HN-1224-TT6**

This product is for research use only and is not intended for diagnostic use.

Product Overview

This kit is performed based on the pre-made Y1H cDNA library, including three major procedures:

1. Bait yeast strain construction;
2. Screen the optimal concentration of AbA for bait strains;
3. Prepare Y1HGold [pBait AbAi] competent state and screen cDNA library.

This product also validates the interaction of positive clones.

Applications

Interaction library screening for nucleic acid

Content

Culture Medium
Self-activation inhibitor
BD plasmid, AD plasmid
Control plasmids
Competent cell

Size

1Set

Storage

Please refer to the manual.