

Product Information

Recombinant Anti-Human MS4A1 Antibody scFv Fragment

Cat. No.: **MOM-18001-S(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody scFv Fragment is against Human CD20, expressed in E. coli

Antigen Description

CD20 is B-lymphocyte surface molecule which plays a role in the development and differentiation of B-cells into plasma cells.

Specific Activity

Tested positive against native antigen.

Target

CD20

Immunogen

Human lymphoblastoid cell line SB.

Source

Mouse

Species Reactivity

Human

Type

scFv Fragment from Mouse IgG1 - kappa

Expression Host

E. coli

Purity

>95%, by SDS-PAGE with silver staining, under reducing conditions.

Applications

Suitable for use in ELISA, WB, Neut and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[MS4A1 membrane-spanning 4-domains, subfamily A, member 1 \[Homo sapiens \]](#)

Official Symbol

MS4A1

Synonyms

MS4A1; membrane-spanning 4-domains, subfamily A, member 1; CD20; B-lymphocyte antigen CD20; B1; Bp35; MS4A2; CD20 antigen; CD20 receptor; leukocyte surface antigen Leu-16; B-lymphocyte cell-surface antigen B1; S7; CVID5; LEU-16; MGC3969;

Gene ID

[931](#)

mRNA Refseq

[NM_021950](#)

Protein Refseq

[NP_068769](#)

MIM

[112210](#)

UniProt ID

P11836

Chromosome Location

11q12-q13.1

Pathway

Hematopoietic cell lineage, organism-specific biosystem; Hematopoietic cell lineage, conserved biosystem;

Function

epidermal growth factor receptor binding;