

Product Information

Recombinant Anti-Human ITGB2 Antibody Fab Fragment

Cat. No.: **MOM-18044-F(E)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody Fab Fragment is against Human ITGB2, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Integrin alpha-L/beta-2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4. Integrins alpha-M/beta-2 and alpha-X/beta-2 are receptors for the iC3b fragment of the third complement component and for fibrinogen. Integrin alpha-X/beta-2 recognizes the sequence G-P-R in fibrinogen alpha-chain. Integrin alpha-M/beta-2 recognizes P1 and P2 peptides of fibrinogen gamma chain. Integrin alpha-M/beta-2 is also a receptor for factor X. Integrin alpha-D/beta-2 is a receptor for ICAM3 and VCAM1.

Specific Activity

Tested positive against native antigen.

Target

ITGB2

Immunogen

The details of the immunogen for this antibody are not available.

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Fab Fragment based on Humanized (from mouse) IgG4

Expression Host

CHO

Predicted N terminal

H chain: QVQLVQS; L Chain: DIQMTQS

Purity

>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords

0

ANTIGEN GENE INFORMATION

Gene Name

[ITGB2 integrin, beta 2 \(complement component 3 receptor 3 and 4 subunit\) \[Homo sapiens \]](#)

Official Symbol

ITGB2

Synonyms

ITGB2; integrin, beta 2 (complement component 3 receptor 3 and 4 subunit); CD18, integrin, beta 2 (antigen CD18 (p95), lymphocyte function associated antigen 1; macrophage antigen 1 (mac 1) beta subunit) , MF17; integrin beta-2; LFA 1; MAC 1; integrin beta chain, beta 2; complement receptor C3 beta-subunit; complement receptor C3 subunit beta; leukocyte cell adhesion molecule CD18; leukocyte-associated antigens CD18/11A, CD18/11B, CD18/11C; cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta; cell surface adhesion glycoprotein LFA-1/CR3/P150,959 beta subunit precursor; LAD; CD18; MF17; MF17; LCAMB; LFA-1; MAC-1;

Gene ID

[3689](#)

mRNA Refseq

[NM_000211](#)

Protein Refseq

[NP_000202](#)

MIM

[600065](#)

UniProt ID

P05107

Chromosome Location

21q22.3

Pathway

Adaptive Immune System, organism-specific biosystem; Amoebiasis, organism-specific biosystem; Amoebiasis, conserved biosystem; CXCR3-mediated signaling events, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem; Cell surface interactions at the vascular wall, organism-specific biosystem;

Function

binding; glycoprotein binding; protein kinase binding; receptor activity;