

Product Information

Recombinant Anti-Human CD22 Antibody Fab Fragment

Cat. No.: **MOM-18047-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody Fab Fragment is bind to Human CD22, expressed in E. coli

Antigen Description

CD22 or cluster of differentiation-22, is a molecule belonging to the SIGLEC family of lectins. It is found on the surface of mature B cells and to a lesser extent on some immature B cells. Generally speaking, CD22 is a regulatory molecule that prevents the overactivation of the immune system and the development of autoimmune diseases. CD22 is a sugar binding transmembrane protein, which specifically binds sialic acid with an immunoglobulin (Ig) domain located at its N-terminus. The presence of Ig domains makes CD22 a member of the immunoglobulin superfamily. CD22 functions as an inhibitory receptor for B cell receptor (BCR) signalling.

Specific Activity

Tested positive against native antigen.

Target

CD22

Immunogen

A combination of recombinant human CD22 extracellular 25 domain and CHO cells engineered to express full-length human CD22 on the cell surface

Source

Mouse

Species Reactivity

Human

Type

Fab Fragment based on Mouse Fab' - G2a - kappa

Expression Host

E. coli

Predicted N terminal

H chain: EVQLVQS; L Chain: EIVLTQS

Purity

>95.0% as determined by Analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[CD22 CD22 molecule \[Homo sapiens \]](#)

Official Symbol

CD22

Synonyms

CD22; CD22 molecule; CD22 antigen; B-cell receptor CD22; sialic acid binding Ig like lectin 2; SIGLEC 2; SIGLEC2; BL-CAM; T-cell surface antigen Leu-14; B-lymphocyte cell adhesion molecule; sialic acid binding Ig-like lectin 2; sialic acid-binding Ig-like lectin 2; SIGLEC-2; FLJ22814; MGC130020;

Gene ID

[933](#)

mRNA Refseq

[NM_001185099](#)

Protein Refseq

[NP_001172028](#)

MIM

[107266](#)

UniProt ID

P20273

Chromosome Location

19q13.1

Pathway

B Cell Receptor Signaling Pathway, organism-specific biosystem; B cell receptor signaling pathway, organism-specific biosystem; B cell receptor signaling pathway, conserved biosystem; BCR signaling pathway, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem; Hematopoietic cell lineage, organism-specific biosystem;

Function

protein binding; sugar binding;