

Product Information

Recombinant Anti-Human IL12+IL23 Antibody Fab Fragment

Cat. No.: **MOM-18063-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Human Antibody Fab Fragment is bind to Human IL12+IL23, expressed in E. coli

Antigen Description

Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Specific Activity

Tested positive against native antigen.

Target

IL12+IL23

Source

Human

Species Reactivity

Human

Type

Fab Fragment based on Human IgG1 - kappa

Expression Host

E. coli

Predicted N terminal

H Chain: EVQLVQS; L Chain: DIQMTQS

Purity

>95.0% as determined by Analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing of samples.

BACKGROUND

Keywords

IL12B; interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40); NKSF2; interleukin-12 subunit beta; CLMF; CLMF2; cytotoxic lymphocyte maturation factor 2; p40; IL 12B; IL12; subunit p40;

interleukin 12; in