

Product Information

Recombinant Anti-Human IGF1R Antibody Fab Fragment

Cat. No.: **MOM-18078-F(E)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Human Antibody Fab Fragment is specific to Human IGF1 Receptor, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Receptor tyrosine kinase which mediates actions of insulin-like growth factor 1 (IGF1). Binds IGF1 with high affinity and IGF2 and insulin (INS) with a lower affinity. The activated IGF1R is involved in cell growth and survival control. IGF1R is crucial for tumor transformation and survival of malignant cell. Ligand binding activates the receptor kinase, leading to receptor autophosphorylation, and tyrosines phosphorylation of multiple substrates, that function as signaling adapter proteins including, the insulin-receptor substrates (IRS1/2), Shc and 14-3-3 proteins.

Specific Activity

Tested positive against native antigen.

Target

IGF1 Receptor

Source

Human

Species Reactivity

Human

Type

Fab Fragment based on Human IgG1 - lambda

Expression Host

CHO

Predicted N terminal

H Chain: EVQLVQS; L Chain: SSELTQD

Purity

>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[IGF1R insulin-like growth factor 1 receptor \[Homo sapiens \]](#)

Official Symbol

IGF1R

Synonyms

IGF1R; insulin-like growth factor 1 receptor; CD221; IGFIR; IGFR; JTK13; MGC18216; IGF-I receptor; soluble IGF1R variant 1; soluble IGF1R variant 2; insulin-like growth factor I receptor; MGC142170; MGC142172;

Gene ID

[3480](#)

mRNA Refseq

[NM_000875](#)

Protein Refseq

[NP_000866](#)

MIM

[147370](#)

UniProt ID

P08069

Chromosome Location

15q26.3

Pathway

Adherens junction, organism-specific biosystem; Adherens junction, conserved biosystem; Apoptosis, organism-specific biosystem; Endochondral Ossification, organism-specific biosystem; Endocytosis, organism-specific biosystem; Endocytosis, conserved biosystem; Focal Adhesion, organism-specific biosystem;

Function

ATP binding; identical protein binding; insulin binding; insulin receptor binding; insulin receptor substrate binding; insulin-like growth factor I binding; insulin-like growth factor binding; insulin-like growth factor-activated receptor activity; nucleotide binding; phosphatidylinositol 3-kinase binding; protein binding; protein tyrosine kinase activity; receptor activity;