

Product Information

Recombinant Anti-Human IL12+IL23 Antibody

Cat. No.: **MOM-18103**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Human Antibody is against Human IL12+IL23, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Specific Activity

Tested positive against native antigen.

Target

IL12+IL23

Immunogen

The details of the immunogen for this antibody are not available.

Source

Human

Species Reactivity

Human

Type

Human IgG1 - lambda

Expression Host

CHO

Predicted N terminal

H chain: QVQLVES; L Chain: QSVLTQP

Purity

>95.0%, determined by analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF, IHC and most other immunological methods.

Storage

4°C. For long term storage, aliquot and store at -20°C. Repeated thawing and freezing must be avoided.

BACKGROUND

Keywords

IL12B; interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40); NKSF2; interleukin-12 subunit beta; CLMF; CLMF2; cytotoxic lymphocyte maturation factor 2; p40; IL 12B; IL12; subunit p40; interleukin 12; in