

Product Information

Recombinant Anti-Human IL12+IL23 Antibody Fab Fragment

Cat. No.: **MOM-18103-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Human Antibody Fab Fragment binds selectively to Human IL12+IL23, expressed in E. coli

Antigen Description

Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.

Specific Activity

Tested positive against native antigen.

Target

IL12+IL23

Immunogen

The details of the immunogen for this antibody are not available.

Source

Human

Species Reactivity

Human

Type

Fab Fragment based on Human IgG1 - lambda

Expression Host

E. coli

Predicted N terminal

H chain: QVQLVES; L Chain: QSVLTQP

Purity

Purity >95% by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

At -20°C for one year.

BACKGROUND

Keywords

IL12B; interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40); NKSF2; interleukin-12 subunit beta; CLMF; CLMF2; cytotoxic lymphocyte maturation factor 2; p40; IL 12B; IL12; subunit p40; interleukin 12; in