

Product Information

Recombinant Anti-Human HLA-DRB1 Antibody Fab Fragment

Cat. No.: **MOM-18150-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody Fab Fragment is specific to Human HLA DR beta, expressed in E. coli

Antigen Description

Human Leukocyte Antigens are highly polymorphic proteins that are involved in the presentation of antigens to the T-cell receptor. There are two classes of HLA antigens, class I (HLA-A, HLA-B and HLA-C) and class II (HLA-D).

Specific Activity

Tested positive against native antigen.

Target

HLA DR beta

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Fab Fragment based on Humanized (from mouse) IgG1

Expression Host

E. coli

Predicted N terminal

H chain: QVQLQES; L Chain: DIQMTQS

Purity

>95.0%, determined by analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[HLA-DRB1 major histocompatibility complex, class II, DR beta 1 \[Homo sapiens \]](#)

Official Symbol

HLA-DRB1

Synonyms

HLA-DRB1; major histocompatibility complex, class II, DR beta 1; HLA DR1B; DW2.2/DR2.2; MHC class II antigen; lymphocyte antigen DRB1; MHC class II HLA-DRw10-beta; human leucocyte antigen DRB1; MHC class II HLA-DR beta 1 chain; MHC class II HLA-DR-beta cell surface glycoprotein; HLA class II histocompatibility antigen, DR-1 beta chain; SS1; DRB1; DRw10; HLA-DRB; HLA-DR1B; FLJ75017; FLJ76359;

Gene ID

[3123](#)

mRNA Refseq

[NM_001243965](#)

Protein Refseq

[NP_001230894](#)

MIM

[142857](#)

UniProt ID

P01911

Chromosome Location

6p21.3

Pathway

Adaptive Immune System, organism-specific biosystem; Allograft rejection, organism-specific biosystem; Allograft rejection, conserved biosystem; Antigen processing and presentation, organism-specific biosystem; Antigen processing and presentation, conserved biosystem; Asthma, organism-specific biosystem; Asthma, conserved biosystem;

Function

MHC class II receptor activity; MHC class II receptor activity;