

Product Information

Recombinant Anti-Human ITGB2 Antibody Fab Fragment

Cat. No.: **MOM-18154-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody Fab Fragment is specific to Human CD18, expressed in E. coli

Antigen Description

Integrin alpha-L/beta-2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4. Integrins alpha-M/beta-2 and alpha-X/beta-2 are receptors for the iC3b fragment of the third complement component and for fibrinogen. Integrin alpha-X/beta-2 recognizes the sequence G-P-R in fibrinogen alpha-chain. Integrin alpha-M/beta-2 recognizes P1 and P2 peptides of fibrinogen gamma chain. Integrin alpha-M/beta-2 is also a receptor for factor X. Integrin alpha-D/beta-2 is a receptor for ICAM3 and VCAM1.

Specific Activity

Tested positive against native antigen.

Target

CD18

Immunogen

Tissue / cell preparation (Human). (Leucocytes of a patient suffering from a LGL-type leukaemia).

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Fab Fragment based on Humanized (from mouse) IgG1

Expression Host

E. coli

Purity

>95%, by SDS-PAGE with silver staining, under reducing conditions.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[ITGB2 integrin, beta 2 \(complement component 3 receptor 3 and 4 subunit\) \[Homo sapiens \]](#)

Official Symbol

ITGB2

Synonyms

ITGB2; integrin, beta 2 (complement component 3 receptor 3 and 4 subunit); CD18, integrin, beta 2 (antigen CD18 (p95), lymphocyte function associated antigen 1; macrophage antigen 1 (mac 1) beta subunit) , MF17; integrin beta-2; LFA 1; MAC 1; integrin beta chain, beta 2; complement receptor C3 beta-subunit; complement receptor C3 subunit beta; leukocyte cell adhesion molecule CD18; leukocyte-associated antigens CD18/11A, CD18/11B, CD18/11C; cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta; cell surface adhesion glycoprotein LFA-1/CR3/P150,959 beta subunit precursor; LAD; CD18; MF17; MF17; LCAMB; LFA-1; MAC-1;

Gene ID

[3689](#)

mRNA Refseq

[NM_000211](#)

Protein Refseq

[NP_000202](#)

MIM

[600065](#)

UniProt ID

P05107

Chromosome Location

21q22.3

Pathway

Adaptive Immune System, organism-specific biosystem; Amoebiasis, organism-specific biosystem; Amoebiasis, conserved biosystem; CXCR3-mediated signaling events, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem; Cell surface interactions at the vascular wall, organism-specific biosystem;

Function

binding; glycoprotein binding; protein kinase binding; receptor activity;