

Product Information

Recombinant Anti-Human IL1B Antibody

Cat. No.: **MOM-18161**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody is specific to Human IL1 beta, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Produced by activated macrophages, IL-1 stimulates thymocyte proliferation by inducing IL-2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammatory response, being identified as endogenous pyrogens, and are reported to stimulate the release of prostaglandin and collagenase from synovial cells.

Specific Activity

Tested positive against native antigen.

Target

IL1 beta

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Humanized (from mouse) IgG2 - kappa

Expression Host

CHO

Predicted N terminal

H chain: QVQLQES; L Chain: DIQMTQS

Purity

>95.0% as determined by analysis by RP-HPLC.

Applications

Suitable for use in IF, IP, Neut, FuncS, ELISA, FC, ICC and most other immunological methods.

Storage

Store at -20°C. Avoid multiple freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[IL1B interleukin 1, beta \[Homo sapiens \]](#)

Official Symbol

IL1B

Synonyms

IL1B; interleukin 1, beta; interleukin-1 beta; IL 1B; IL1 BETA; IL1F2; IL-1 beta; catabolin; preinterleukin 1 beta; pro-interleukin-1-beta; IL-1; IL1-BETA;

Gene ID

[3553](#)

mRNA Refseq

[NM_000576](#)

Protein Refseq

[NP_000567](#)

MIM

[147720](#)

UniProt ID

P01584

Chromosome Location

2q14

Pathway

African trypanosomiasis, organism-specific biosystem; African trypanosomiasis, conserved biosystem; Alzheimers disease, organism-specific biosystem; Alzheimers disease, conserved biosystem; Amoebiasis, organism-specific biosystem; Amoebiasis, conserved biosystem; Apoptosis, organism-specific biosystem;

Function

cytokine activity; cytokine activity; growth factor activity; interleukin-1 receptor binding; protein domain specific binding;