

Product Information

Recombinant Anti-Human CD4 Antibody

Cat. No.: **MOM-18171**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Chimeric (primate/human) Antibody binds selectively to Human CD4, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Accessory protein for MHC class-II antigen/T-cell receptor interaction. May regulate T-cell activation. Induces the aggregation of lipid rafts.

Specific Activity

Tested positive against native antigen.

Target

CD4

Immunogen

Recombinant fragment (Human) corresponding to the external domain.

Source

Chimeric (primate/human)

Species Reactivity

Human

Type

Chimeric (primate/human) IgG1 - lambda

Expression Host

CHO

Purity

>95.0% as determined by analysis by SDS-PAGE.

Applications

Suitable for use in IP, IF, FuncS, FC, Neut, ELISA and most other immunological methods.

Storage

Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[CD4 CD4 molecule \[Homo sapiens \]](#)

Official Symbol

CD4

Synonyms

CD4; CD4 molecule; CD4 antigen (p55) , T cell surface glycoprotein CD4; T-cell surface glycoprotein CD4; CD4 receptor; CD4 antigen (p55); T-cell surface antigen T4/Leu-3; CD4mut;

Gene ID

[920](#)

mRNA Refseq

[NM_000616](#)

Protein Refseq

[NP_000607](#)

MIM

[186940](#)

UniProt ID

P01730

Chromosome Location

12p13.31

Pathway

Adaptive Immune System, organism-specific biosystem; Alpha-defensins, organism-specific biosystem; Antigen processing and presentation, organism-specific biosystem; Antigen processing and presentation, conserved biosystem; Arf1 pathway, organism-specific biosystem; Binding and entry of HIV virion, organism-specific biosystem; C-MYB transcription factor network, organism-specific biosystem;

Function

MHC class II protein binding; coreceptor activity; enzyme binding; extracellular matrix structural constituent; glycoprotein binding; protein binding; protein homodimerization activity; protein kinase binding; receptor activity; transmembrane signaling receptor activity; zinc ion binding;