

Product Information

Recombinant Anti-Human CCR4 Antibody Fab Fragment

Cat. No.: **MOM-18187-F(E)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody Fab Fragment is bind to Human CCR4, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

High affinity receptor for the C-C type chemokines CCL17/TARC and CCL22/MDC. The activity of this receptor is mediated by G(i) proteins which activate a phosphatidylinositol-calcium second messenger system. Can function as a chemoattractant homing receptor on circulating memory lymphocytes and as a coreceptor for some primary HIV-2 isolates. In the CNS, could mediate hippocampal-neuron survival.

Specific Activity

Tested positive against native antigen.

Target

CCR4

Immunogen

KLH conjugate CCR4 partial peptide

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Fab Fragment based on Humanized (from mouse) IgG1 - kappa

Expression Host

CHO

Predicted N terminal

H chain: EVQLVES; L Chain: DIVMTQS

Purity

>95.0%. Determined by analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing of samples.

ANTIGEN GENE INFORMATION

Gene Name

[CCR4 chemokine \(C-C motif\) receptor 4 \[Homo sapiens \]](#)

Official Symbol

CCR4

Synonyms

CCR4; chemokine (C-C motif) receptor 4; C-C chemokine receptor type 4; CC CKR 4; CD194; ChemR13; CKR4; CMKBR4; k5 5; CCR-4; C-C CKR-4; chemokine (C-C) receptor 4; K5-5; CC-CKR-4; HGCN:14099; MGC88293;

Gene ID

[1233](#)

mRNA Refseq

[NM_005508](#)

Protein Refseq

[NP_005499](#)

MIM

[604836](#)

UniProt ID

P51679

Chromosome Location

3p24-p21.3

Pathway

Chemokine receptors bind chemokines, organism-specific biosystem; Chemokine signaling pathway, organism-specific biosystem; Chemokine signaling pathway, conserved biosystem; Class A/1 (Rhodopsin-like receptors), organism-specific biosystem; Cytokine-cytokine receptor interaction, organism-specific biosystem; Cytokine-cytokine receptor interaction, conserved biosystem; G alpha (i) signalling events, organism-specific biosystem;

Function

C-C chemokine receptor activity; G-protein coupled receptor activity; chemokine receptor activity; receptor activity; signal transducer activity;