

Product Information

Recombinant Anti-Human SDC1 Antibody Fab Fragment

Cat. No.: **MOM-18197-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Chimeric (mouse/human) Antibody Fab Fragment is against Human SDC1, expressed in E. coli

Antigen Description

The protein encoded by this gene is a transmembrane (type I) heparan sulfate proteoglycan and is a member of the syndecan proteoglycan family. The syndecans mediate cell binding, cell signaling, and cytoskeletal organization and syndecan receptors are required for internalization of the HIV-1 tat protein. The syndecan-1 protein functions as an integral membrane protein and participates in cell proliferation, cell migration and cell-matrix interactions via its receptor for extracellular matrix proteins. Altered syndecan-1 expression has been detected in several different tumor types. While several transcript variants may exist for this gene, the full-length nature of only two have been described to date. These two represent the major variants of this gene and encode the same protein.

Specific Activity

Tested positive against native antigen.

Target

SDC1

Immunogen

U266 cell line

Source

Chimeric (mouse/human)

Species Reactivity

Human

Type

Fab Fragment based on Chimeric (mouse/human) IgG4 - kappa

Expression Host

E. coli

Purity

>95.0% as determined by analysis by RP-HPLC.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[SDC1 syndecan 1 \[Homo sapiens\]](#)

Official Symbol

SDC1

Synonyms

SDC1; syndecan 1; SDC; syndecan-1; CD138; SYND1; syndecan; syndecan proteoglycan 1; CD138 antigen; heparan sulfate proteoglycan fibroblast growth factor receptor;

Gene ID

[6382](#)

mRNA Refseq

[NM_001006946](#)

Protein Refseq

[NP_001006947](#)

MIM

[186355](#)

UniProt ID

P18827

Chromosome Location

2p24.1

Pathway

Arf6 trafficking events, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem; Chylomicron-mediated lipid transport, organism-specific biosystem; ECM-receptor interaction, organism-specific biosystem; ECM-receptor interaction, conserved biosystem; Integrins in angiogenesis, organism-specific biosystem;

Function

cytoskeletal protein binding; glycoprotein binding; lipoprotein lipase activity; protein C-terminus binding; protein binding;