

Product Information

Recombinant Anti-Human MSLN Antibody scFv Fragment

Cat. No.: **MOM-18201-S(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Chimeric (mouse/human) Antibody scFv Fragment is directed against Human Mesothelin, expressed in *E. coli*

Antigen Description

Membrane-anchored forms may play a role in cellular adhesion. Megakaryocyte-potentiating factor (MPF) potentiates megakaryocyte colony formation in vitro.

Specific Activity

Tested positive against native antigen.

Target

Mesothelin

Immunogen

Spleen lymphocytes from mice fused with OVCAR-3, an ovarian tumor cell line derived from poorly differentiated adenocarcinoma of the ovary.

Source

Chimeric (mouse/human)

Species Reactivity

Human

Type

scFv Fragment from Chimeric (mouse/human) IgG1 - kappa

Expression Host

E. coli

Purity

>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.

Applications

Suitable for use in ELISA, WB, Neut and most other immunological methods.

Storage

Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

ANTIGEN GENE INFORMATION

Gene Name

[MSLN mesothelin \[Homo sapiens \]](#)

Official Symbol

MSLN

Synonyms

MSLN; mesothelin; CAK1; MPF; CAK1 antigen; megakaryocyte potentiating factor; soluble MPF mesothelin related protein; pre-pro-megakaryocyte-potentiating factor; SMRP;

Gene ID

[10232](#)

mRNA Refseq

[NM_001177355](#)

Protein Refseq

[NP_001170826](#)

MIM

[601051](#)

UniProt ID

Q13421

Chromosome Location

16p13.3