

Product Information

Recombinant Anti-Human alcam Antibody Fab Fragment

Cat. No.: **MOM-18268-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody Fab Fragment binds selectively to Human ALCAM, expressed in E. coli

Antigen Description

Cell adhesion molecule that binds to CD6. Involved in neurite extension by neurons via heterophilic and homophilic interactions. May play a role in the binding of T- and B-cells to activated leukocytes, as well as in interactions between cells of the nervous system.

Specific Activity

Tested positive against native antigen.

Target

ALCAM

Immunogen

Purified recombinant fragment of Human CD166 (aa405-524) expressed in E. Coli.

Source

Mouse

Species Reactivity

Human

Type

Fab

Expression Host

E. coli

Purity

Purity >95% by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[ALCAM activated leukocyte cell adhesion molecule \[Homo sapiens \]](#)

Official Symbol

ALCAM

Synonyms

ALCAM; activated leukocyte cell adhesion molecule; activated leucocyte cell adhesion molecule; CD166 antigen; CD166; MEMD; FLJ38514; MGC71733

Gene ID

[214](#)

mRNA Refseq

[NM_001243280](#)

Protein Refseq

[NP_001230209](#)

MIM

[601662](#)

UniProt ID

Q13740

Chromosome Location

3q13.1

Pathway

Axon guidance, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem; Developmental Biology, organism-specific biosystem; L1CAM interactions, organism-specific biosystem;

Function

receptor binding;