

Product Information

Recombinant Anti-Human acvrl1 Antibody

Cat. No.: **MOM-18270**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody binds selectively to Human ACVRL1, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators. Receptor for TGF-beta. May bind activin as well.

Specific Activity

Tested positive against native antigen.

Target

ACVRL1

Source

Mouse

Species Reactivity

Human

Type

IgG

Expression Host

CHO

Purity

>95.0% as determined by analysis by SDS-PAGE.

Applications

Suitable for use in IF, IP, Neut, FuncS, ELISA, FC, IHC and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[ACVRL1 activin A receptor type II-like 1 \[Homo sapiens \]](#)

Official Symbol

ACVRL1

Synonyms

ACVRL1; activin A receptor type II-like 1; ACVRLK1, ORW2; serine/threonine-protein kinase receptor R3; ALK1; HHT; HHT2; activin receptor-like kinase 1; TGF-B superfamily receptor type I; activin A receptor, type II-like kinase 1; ORW2; SKR3; ALK-1; TSR-I; ACVRLK1

Gene ID

[94](#)

mRNA Refseq

[NM_000020](#)

Protein Refseq

[NP_000011](#)

MIM

[601284](#)

UniProt ID

P37023

Chromosome Location

12q11-q14

Pathway

ALK1 pathway, organism-specific biosystem; ALK1 signaling events, organism-specific biosystem; Id Signaling Pathway, organism-specific biosystem; TGF-beta Receptor Signaling Pathway, organism-specific biosystem;

Function

ATP binding; SMAD binding; activin binding; activin receptor activity, type I; contributes_to activin receptor activity, type I; metal ion binding; nucleotide binding; protein binding; protein serine/threonine kinase activity; receptor activity; receptor signaling protein serine/threonine kinase activity; transforming growth factor beta binding; contributes_to transforming growth factor beta receptor activity, type I; transforming growth factor beta-activated receptor activity; transmembrane receptor protein serine/threonine kinase activity;