

# Product Information

## Recombinant Anti-Human lypd3 Antibody

Cat. No.: **MOM-18296**

This product is for research use only and is not intended for diagnostic use.

### Product Overview

Recombinant Mouse Antibody is against Human LYPD3, expressed in Chinese Hamster Ovary cells(CHO)

### Antigen Description

Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.

### Specific Activity

Tested positive against native antigen.

### Target

LYPD3

### Source

Mouse

### Species Reactivity

Human

### Type

IgG

### Expression Host

CHO

### Purity

>95%, by SDS-PAGE with silver staining, under reducing conditions.

### Applications

Suitable for use in Neut and most other immunological methods.

### Storage

Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.

## ANTIGEN GENE INFORMATION

### Gene Name

[LYPD3 LY6/PLAUR domain containing 3 \[ Homo sapiens \]](#)

### Official Symbol

LYPD3

### Synonyms

LYPD3; LY6/PLAUR domain containing 3; ly6/PLAUR domain-containing protein 3; C4.4A; MIG-C4; 2310061G07Rik; matrigel-induced gene C4 protein; GPI-anchored metastasis-associated protein homolog; GPI-anchored metastasis-associated protein C4.4A homolog

**Gene ID**

[27076](#)

**mRNA Refseq**

[NM\\_014400](#)

**Protein Refseq**

[NP\\_055215](#)

**MIM**

[609484](#)

**UniProt ID**

O95274

**Chromosome Location**

19q13.31