

Product Information

Recombinant Anti-Human lypd3 Antibody Fab Fragment

Cat. No.: **MOM-18296-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody Fab Fragment binds selectively to Human LYPD3, expressed in E. coli

Antigen Description

Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.

Specific Activity

Tested positive against native antigen.

Target

LYPD3

Source

Mouse

Species Reactivity

Human

Type

Fab

Expression Host

E. coli

Purity

>95.0% as determined by analysis by RP-HPLC.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C. Avoid multiple freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[LYPD3 LY6/PLAUR domain containing 3 \[Homo sapiens \]](#)

Official Symbol

LYPD3

Synonyms

LYPD3; LY6/PLAUR domain containing 3; ly6/PLAUR domain-containing protein 3; C4.4A; MIG-C4; 2310061G07Rik; matrigel-induced gene C4 protein; GPI-anchored metastasis-associated protein homolog; GPI-anchored metastasis-associated protein C4.4A homolog

Gene ID

[27076](#)

mRNA Refseq

[NM_014400](#)

Protein Refseq

[NP_055215](#)

MIM

[609484](#)

UniProt ID

O95274

Chromosome Location

19q13.31