

Product Information

Recombinant Anti-Human met Antibody

Cat. No.: **MOM-18348**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody is specific to Human MET, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Receptor for hepatocyte growth factor and scatter factor. Has a tyrosine-protein kinase activity. Functions in cell proliferation, scattering, morphogenesis and survival.

Specific Activity

Tested positive against native antigen.

Target

MET

Immunogen

Synthetic peptide (KLH coupled) corresponding to residues surrounding Tyr1234 of human Met.

Source

Mouse

Species Reactivity

Human

Type

IgG

Expression Host

CHO

Purity

>95.0% as determined by Analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut and most other immunological methods.

Storage

Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[MET met proto-oncogene \(hepatocyte growth factor receptor\) \[Homo sapiens \]](#)

Official Symbol

MET

Synonyms

MET; met proto-oncogene (hepatocyte growth factor receptor); hepatocyte growth factor receptor; HGFR; RCCP2; SF receptor; HGF receptor; oncogene MET; HGF/SF receptor; proto-oncogene c-Met; scatter factor receptor; tyrosine-protein kinase Met; met proto-oncogene tyrosine kinase; AUTS9; c-Met

Gene ID

[4233](#)

mRNA Refseq

[NM_000245](#)

Protein Refseq

[NP_000236](#)

MIM

[164860](#)

UniProt ID

P08581

Chromosome Location

7q31

Pathway

Adherens junction, organism-specific biosystem; Adherens junction, conserved biosystem; Alpha6-Beta4 Integrin Signaling Pathway, organism-specific biosystem; Arf6 signaling events, organism-specific biosystem; Axon guidance, organism-specific biosystem; Axon guidance, conserved biosystem; Axon guidance, organism-specific biosystem;

Function

ATP binding; hepatocyte growth factor-activated receptor activity; nucleotide binding; protein binding; protein tyrosine kinase activity; protein tyrosine kinase activity; receptor activity;