

Product Information

Recombinant Anti-Human tac3 Antibody Fab Fragment

Cat. No.: **MOM-18446-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody Fab Fragment is specific to Human TAC3, expressed in E. coli

Antigen Description

Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodilators and secretagogues, and contract (directly or indirectly) many smooth muscles.

Specific Activity

Tested positive against native antigen.

Target

TAC3

Source

Mouse

Species Reactivity

Human

Type

Fab

Expression Host

E. coli

Purity

>95.0% as determined by analysis by RP-HPLC.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C. Avoid multiple freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[TAC3 tachykinin 3 \[Homo sapiens \]](#)

Official Symbol

TAC3

Synonyms

TAC3; tachykinin 3; neurokinin beta , neuromedin K , NKNB; tachykinin-3; NKB; ZNEUROK1; neuromedin K; gamma tachykinin 3; preprotachykinin B; NKNB; PRO1155

Gene ID

[6866](#)

mRNA Refseq

[NM_001178054](#)

Protein Refseq

[NP_001171525](#)

MIM

[162330](#)

UniProt ID

Q9UHF0

Chromosome Location

12q13-q21

Pathway

Class A/1 (Rhodopsin-like receptors), organism-specific biosystem; G alpha (q) signalling events, organism-specific biosystem; GPCR downstream signaling, organism-specific biosystem; GPCR ligand binding, organism-specific biosystem; Peptide ligand-binding receptors, organism-specific biosystem; Signal Transduction, organism-specific biosystem; Signaling by GPCR, organism-specific biosystem;

Function

receptor binding;