

Product Information

Recombinant Anti-Human *lilrb3* Antibody

Cat. No.: **MOM-18471**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody is against Human LILRB3, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

LILRB3 is a 631 amino acid type I transmembrane glycoprotein, which contains four immunoreceptor tyrosine-based inhibition motif (ITIM) sequences within a long cytoplasmic tail. Phosphorylation of the tyrosine residues within ITIMs is known to enable the binding and activation of protein tyrosine phosphatases, which act as cell signalling modulators and inhibitors of cell activation. LILRB3 may act as receptor for class I MHC antigens.

Specific Activity

Tested positive against native antigen.

Target

LILRB3

Source

Mouse

Species Reactivity

Human

Type

IgG

Expression Host

CHO

Purity

>95.0%. Determined by analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in Neut, ICC and most other immunological methods.

Storage

Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.

ANTIGEN GENE INFORMATION

Gene Name

[LILRB3 leukocyte immunoglobulin-like receptor, subfamily B \(with TM and ITIM domains\), member 3 \[Homo sapiens \]](#)

Official Symbol

LILRB3

Synonyms

LILRB3; leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 3; leukocyte immunoglobulin-like receptor subfamily B member 3; CD85a; HL9; ILT5; LIR 3; LIR3; immunoglobulin-like transcript 5; monocyte inhibitory receptor HL9; CD85 antigen-like family member A; leukocyte immunoglobulin-like receptor 3; PIRB; CD85A; ILT-5; LIR-3; MGC138403

Gene ID

[11025](#)

mRNA Refseq

[NM_001081450](#)

Protein Refseq

[NP_001074919](#)

MIM

[604820](#)

UniProt ID

O75022

Chromosome Location

19q13.4

Pathway

Adaptive Immune System, organism-specific biosystem; B cell receptor signaling pathway, organism-specific biosystem; B cell receptor signaling pathway, conserved biosystem; Immune System, organism-specific biosystem; Immunoregulatory interactions between a Lymphoid and a non-Lymphoid cell, organism-specific biosystem; Osteoclast differentiation, organism-specific biosystem; Osteoclast differentiation, conserved biosystem;

Function

protein binding; receptor activity; transmembrane signaling receptor activity;