

Product Information

Recombinant Anti-Human gpa33 Antibody

Cat. No.: **MOM-18569**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody is against Human GPA33, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

May play a role in cell-cell recognition and signaling.

Specific Activity

Tested positive against native antigen.

Target

GPA33

Immunogen

The details of the immunogen for this antibody are not available.

Source

Mouse

Species Reactivity

Human

Type

IgG

Expression Host

CHO

Purity

>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.

Applications

Suitable for use in Neut, IF, ICC and most other immunological methods.

Storage

Store at -20°C. Avoid multiple freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[GPA33 glycoprotein A33 \(transmembrane\) \[Homo sapiens \]](#)

Official Symbol

GPA33

Synonyms

GPA33; glycoprotein A33 (transmembrane); cell surface A33 antigen; A33; transmembrane glycoprotein A33; MGC129986; MGC129987;

Gene ID

[10223](#)

mRNA Refseq

[NM_005814](#)

Protein Refseq

[NP_005805](#)

MIM

[602171](#)

UniProt ID

Q99795

Chromosome Location

1q24.1

Function

receptor activity;