

Product Information

Recombinant Anti-Human itgav Antibody scFv Fragment

Cat. No.: **MOM-18583-S(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody scFv Fragment is directed against Human ITGAV, expressed in E. coli

Antigen Description

The alpha-V integrins are receptors for vitronectin, cytotactin, fibronectin, fibrinogen, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin and vWF. They recognize the sequence R-G-D in a wide array of ligands. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions.

Specific Activity

Tested positive against native antigen.

Target

ITGAV

Immunogen

Full length native protein (purified) (Human).

Source

Mouse

Species Reactivity

Human

Type

scFv

Expression Host

E. coli

Purity

>95.0% as determined by Analysis by RP-HPLC & analysis by SDS-PAGE.

Applications

Suitable for use in ELISA, WB, Neut and most other immunological methods.

Storage

Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

ANTIGEN GENE INFORMATION

Gene Name

[ITGAV integrin, alpha V \[Homo sapiens \]](#)

Official Symbol

ITGAV

Synonyms

ITGAV; integrin, alpha V; antigen identified by monoclonal L230 , integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51) , MSK8, vitronectin receptor , VNRA, VTNR; integrin alpha-V; CD51; integrin alphaVbeta3; vitronectin receptor subunit alpha; antigen identified by monoclonal L230; integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51); MSK8; VNRA; VTNR; DKFZp686A08142;

Gene ID

[3685](#)

mRNA Refseq

[NM_001144999](#)

Protein Refseq

[NP_001138471](#)

MIM

[193210](#)

UniProt ID

P06756

Chromosome Location

2q31-q32

Pathway

Adaptive Immune System, organism-specific biosystem; Antigen processing-Cross presentation, organism-specific biosystem; Arrhythmogenic right ventricular cardiomyopathy (ARVC), organism-specific biosystem; Arrhythmogenic right ventricular cardiomyopathy (ARVC), conserved biosystem; Axon guidance, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem;

Function

contributes_to insulin-like growth factor I binding; contributes_to opsonin binding; protein binding; protein kinase C alpha binding; contributes_to protein kinase C binding; receptor activity; transforming growth factor beta binding;