

Product Information

Recombinant Anti-Human tnc Antibody Fab Fragment

Cat. No.: **MOM-18623-F(E)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Mouse Antibody Fab Fragment is bind to Human TNC, expressed in Chinese Hamster Ovary cells(CHO)

Antigen Description

Extracellular matrix protein implicated in guidance of migrating neurons as well as axons during development, synaptic plasticity as well as neuronal regeneration. Promotes neurite outgrowth from cortical neurons grown on a monolayer of astrocytes. Ligand for integrins alpha-8/beta-1, alpha-9/beta-1, alpha-V/beta-3 and alpha-V/beta-6.

Specific Activity

Tested positive against native antigen.

Target

TNC

Immunogen

Full length native protein.

Source

Mouse

Species Reactivity

Human

Type

Fab

Expression Host

CHO

Purity

>95.0% as determined by analysis by SDS-PAGE.

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[TNC tenascin C \[Homo sapiens \]](#)

Official Symbol

TNC

Synonyms

TNC; tenascin C; hexabrachion (tenascin C, cytotactin) , HXB; tenascin; hexabrachion (tenascin); MGC167029; TN; GP 150-225; cytotactin; neuronectin; myotendinous antigen; tenascin-C isoform 14/AD1/16; glioma-associated-extracellular matrix antigen; GP; JI; HXB; GMEM; TN-C; 150-225;

Gene ID

[3371](#)

mRNA Refseq

[NM_002160](#)

Protein Refseq

[NP_002151](#)

MIM

[187380](#)

UniProt ID

P24821

Chromosome Location

9q32-q34

Pathway

ECM-receptor interaction, organism-specific biosystem; ECM-receptor interaction, conserved biosystem; Focal Adhesion, organism-specific biosystem; Focal adhesion, organism-specific biosystem; Focal adhesion, conserved biosystem; Integrin cell surface interactions, organism-specific biosystem; Signal Transduction, organism-specific biosystem;

Function

binding; receptor binding; syndecan binding;