

Product Information

Recombinant Anti-Human ITGAV Antibody Fab Fragment

Cat. No.: **MOM-H01-F(P)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody Fab Fragment is directed against Human ITGAV, expressed in E. coli

Antigen Description

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

Specific Activity

ITGAV (integrin alpha V, CD51) [Homo sapiens]

Target

ITGAV

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Humanized (from mouse) Fab-IgG2 - kappa

Expression Host

E. coli

Purity

>95.0% as determined by analysis by RP-HPLC.

Purification

Purified by Nickel ion affinity chromatography

Applications

Suitable for use in FC, IP, ELISA, Neut, FuncS, IF and most other immunological methods.

Cellular Localization

kappa

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month. Avoid freeze/thaw cycles.

ANTIGEN GENE INFORMATION

Gene Name

[ITGAV integrin, alpha V \[Homo sapiens \]](#)

Official Symbol

ITGAV

Synonyms

ITGAV; integrin, alpha V; antigen identified by monoclonal L230 , integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51) , MSK8, vitronectin receptor , VNRA, VTNR; integrin alpha-V; CD51; integrin alphaVbeta3; vitronectin receptor subunit alpha; antigen identified by monoclonal L230; integrin, alpha V (vitronectin receptor, alpha polypeptide, antigen CD51); MSK8; VNRA; VTNR; DKFZp686A08142;

Gene ID

[3685](#)

mRNA Refseq

[NM_001144999](#)

Protein Refseq

[NP_001138471](#)

MIM

[193210](#)

UniProt ID

P06756

Chromosome Location

2q31-q32

Pathway

Adaptive Immune System, organism-specific biosystem; Antigen processing-Cross presentation, organism-specific biosystem; Arrhythmogenic right ventricular cardiomyopathy (ARVC), organism-specific biosystem; Arrhythmogenic right ventricular cardiomyopathy (ARVC), conserved biosystem; Axon guidance, organism-specific biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem;

Function

contributes_to insulin-like growth factor I binding; contributes_to opsonin binding; protein binding; protein kinase C alpha binding; contributes_to protein kinase C binding; receptor activity; transforming growth factor beta binding;