

Product Information

Recombinant Anti-Human IL17A Antibody Fab Fragment

Cat. No.: **MOM-H57-F(E)**

This product is for research use only and is not intended for diagnostic use.

Product Overview

Recombinant Humanized (from mouse) Antibody Fab Fragment is specific for Human IL17A, expressed in HEK293

Antigen Description

The protein encoded by this gene is a proinflammatory cytokine produced by activated T cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and cyclooxygenase-2

Specific Activity

IL17A (interleukin 17A, IL-17A) [Homo sapiens]

Target

IL17A

Source

Humanized (from mouse)

Species Reactivity

Human

Type

Humanized (from mouse) Fab-IgG1 - kappa

Expression Host

HEK293

Purity

>95.0% as determined by analysis by RP-HPLC.

Purification

Purified by Nickel ion affinity chromatography

Applications

Suitable for use in Neut, ELISA and most other immunological methods.

Cellular Localization

kappa

Storage

Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

ANTIGEN GENE INFORMATION

Gene Name

[IL17A interleukin 17A \[Homo sapiens \]](#)

Official Symbol

IL17A

Synonyms

IL17A; interleukin 17A; CTLA8, IL17, interleukin 17 (cytotoxic T lymphocyte associated serine esterase 8); interleukin-17A; cytotoxic T lymphocyte associated protein 8; IL 17; IL 17A; CTLA-8; cytotoxic T-lymphocyte-associated antigen 8; cytotoxic T-lymphocyte-associated protein 8; cytotoxic T-lymphocyte-associated serine esterase 8; interleukin 17 (cytotoxic T-lymphocyte-associated serine esterase 8); IL17; CTLA8; IL-17; IL-17A;

Gene ID

[3605](#)

mRNA Refseq

[NM_002190](#)

Protein Refseq

[NP_002181](#)

MIM

[603149](#)

UniProt ID

Q16552

Chromosome Location

6p12

Pathway

Cytokine-cytokine receptor interaction, organism-specific biosystem; Cytokine-cytokine receptor interaction, conserved biosystem; IL23-mediated signaling events, organism-specific biosystem; IL27-mediated signaling events, organism-specific biosystem; Rheumatoid arthritis, organism-specific biosystem; Rheumatoid arthritis, conserved biosystem;

Function

cytokine activity;