

Product Information

MemDX™ Membrane Protein Human TMEM35 (Transmembrane protein 35A)

Cat. No.: **MP0020J**

This product is for research use only and is not intended for diagnostic use.

This product is a 18.3 kDa Human TMEM35 membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM35

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

18.3 kDa

TMD

4

Sequence

MASPRVTIVALSVLGLFFVFMGTIKLTPRLSKDAYSEMKRAYKSYVRALPLLKKMGINSILLRKSIGA
LEVACGIVMTLVPGRPKDVANFFLLLLVLAVLFFHQLVGDPLKRYAHALVFGILLTCRLLIARKPEDRSS
EKKPLPGNAEEQPSLYEKAPQGKVKVS

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Powder

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

TMEM35

Full Name

Transmembrane protein 35A

Introduction

A soluble peptide released by shedding may interact with NGFR and modulate neurite outgrowth.

Alternative Names

FLJ14084; TMEM35;NACHO; TUF-1; nAChR regulator; nicotinic acetylcholine receptor regulator

Gene ID

[59353](#)

UniProt ID

[Q53FP2](#)