

Product Information

MemDX™ Membrane Protein Human IL12B (Interleukin 12B) expressed in Hi-5 insect for Antibody Discovery

Cat. No.: **MP0034Q**

This product is for research use only and is not intended for diagnostic use.

This product is a 80 kDa Human IL12B membrane protein expressed in Hi-5 insect. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

IL12B

Protein Length

Partial

Protein Class

Druggable Genome, Secreted Protein, Transmembrane

Molecular Weight

80 kDa

Sequence

IWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEVLGSGKTLTIQVKEFGDAGQYTCHKSGGEVLSHSLLLHKKKE

Product Description

Expression Systems

Hi-5 insect

Tag

His

Form

Powder

Endotoxin

< 1 EU/μg

Purity

>95% as determined by SDS-PAGE and Coomassie blue staining

Buffer

0.2 µM filtered solution of 20mM phosphate buffer, 100mM NaCl, pH 7.2

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

IL12B

Full Name

Interleukin 12B

Introduction

This gene encodes a subunit of interleukin 12, a cytokine that acts on T and natural killer cells, and has a broad array of biological activities. Interleukin 12 is a disulfide-linked heterodimer composed of the 40 kD cytokine receptor like subunit encoded by this gene, and a 35 kD subunit encoded by IL12A. This cytokine is expressed by activated macrophages that serve as an essential inducer of Th1 cells development. This cytokine has been found to be important for sustaining a sufficient number of memory/effector Th1 cells to mediate long-term protection to an intracellular pathogen.

Overexpression of this gene was observed in the central nervous system of patients with multiple sclerosis (MS), suggesting a role of this cytokine in the pathogenesis of the disease. The promoter polymorphism of this gene has been reported to be associated with the severity of atopic and non-atopic asthma in children.

Alternative Names

CLMF; CLMF2; IL-12B; IMD28; IMD29; NKSF; NKSF2; interleukin-12 subunit beta; CLMF p40; NK cell stimulatory factor chain 2; IL12, subunit p40; cytotoxic lymphocyte maturation factor 40 kDa subunit; interleukin 12, p40; interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40); interleukin-12 beta chain; natural killer cell stimulatory factor, 40 kD subunit

Gene ID

[3593](#)

UniProt ID

[P29460](#)