

Product Information

MemDX™ Membrane Protein Human APCDD1 (APC down-regulated 1) with His-Flag tag for Antibody Discovery

Cat. No.: **MP0044X**

This product is for research use only and is not intended for diagnostic use.

This product is a 57.5 kDa Human APCDD1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

APCDD1

Molecular Weight

57.5 kDa

Product Description

Application

SDS-PAGE

Expression Systems

HEK293

Tag

His-Flag

Form

Liquid

Endotoxin

< 1.0 EU/ug

Purification

Ni affinity column

Purity

> 90% as determined by SEC-HPLC

Buffer

20 mM Phosphate, 150 mM NaCl, pH7.8 (10% Glycerol)

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

APCDD1

Full Name

APC down-regulated 1

Introduction

This locus encodes an inhibitor of the Wnt signaling pathway. Mutations at this locus have been associated with hereditary hypotrichosis simplex. Increased expression of this gene may also be associated with colorectal carcinogenesis

Alternative Names

B7323; DRAPC1; FP7019

Gene ID

[147495](#)

UniProt ID

[Q8J025](#)