

Product Information

MemDX™ Membrane Protein Human TMEM184C (Transmembrane protein 184C)

Cat. No.: **MP0045J**

This product is for research use only and is not intended for diagnostic use.

This product is a 50 kDa Human TMEM184C membrane protein expressed in HEK293T. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM184C

Protein Length

Full-length

Protein Class

Transmembrane

Molecular Weight

50 kDa

TMD

7

Sequence

MPCTCTWRNWRQWIRPLVAVIYLVSIIVAVPLCVWELQKLEVGIHTKAWFIAGIFLLLTIPISLWVILQH
LVHYTQPELQKPIIRILWMVPIYSLDSWIALKYPGIAIYVDTCRECYEAYVIYNFMGFLTNYLTNRYPNL
VLILEAKDQQKHFPPLCCPPWAMGEVLLFRCKLGVLYQYTVVRPFTTIVALICELLGIYDEGNFSFSNAW
TYLVIINNMSQLFAMYCLLLFYKVLKEELSPIQPVGKFLCVKLVVFSFWQAVVIALLVKVGVISEKHTW
EWQTVEAVATGLQDFIICIEMLAAIAHHYTF SYKPYVQEAEEGSCFDSFLAMWDVSDIRDDISEQVRHV
GRTVRGHPRKLPEDQDQNEHTSLLSSSSQDAISIASSMPPSPMGHYQGFGHTVTPQTPTTAKISDEI
LSDTIGEKKEPSDKSVDS

Product Description

Expression Systems

HEK293T

Tag

C-Myc/DDK

Form

Liquid

Purification

Anti-DDK affinity column followed by conventional chromatography steps

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer

25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at +4°C for up to one week or several months at -80°C

Target**Target Protein**

TMEM184C

Full Name

Transmembrane protein 184C

Introduction

Possible tumor suppressor which may play a role in cell growth.

Alternative Names

TMEM34; transmembrane protein 34

Gene ID

[55751](#)

UniProt ID

[Q9NVA4](#)